


Development of an SPE-LC/MS/MS Assay for the Simultaneous Quantification of Amyloid Beta Peptides in Cerebrospinal Fluid in Support of Alzheimer's Research

Dr. Erin E. Chambers
Waters Corporation

Presented by Dr. Diego Rodriguez Cabaleiro
Waters Europe

Background: Amyloid β Peptides

- Clinical significance of amyloid β peptides
 - Formation of aggregates/plaques in the brain is a critical event in Alzheimer's Disease
 - Drug development strategies aimed at lowering production of these peptides or enhancing their clearance
 - Inhibition/modulation of beta and gamma secretase enzyme
- Current analytical methodologies
 - ELISA- based assays
 - Immunoprecipitation (IP) followed by LCMSMS



Amyloid β 1-38

DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGG

MW 4132, pI 5.2, HPLC index 96

Amyloid β 1-40

DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGV

MW 4330, pI 5.2, HPLC index 103

Amyloid β 1-42

DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVIA

MW 4516, pI 5.2, HPLC index 117

Internal standards are N15 labeled versions of 1-38, 1-40, and 1-42

Waters Regulated Bioanalysis System Solution

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Sample Preparation Solutions

- Best in class

ACQUITY UPLC® I-Class

- The fastest and most resolving LC without compromise in robustness

ACQUITY UPLC Columns

- Different choices for best performance

® TQ-S

- Most sensitive MS

UNIFI

- Compliant, Interactive workflow-driven data platform

Service


- Installation, maintenance & training
- Compliance services



- Introduction
- Mass Spectrometry
- Liquid Chromatography
- Choice of Standard Curve Matrix
- Solid Phase Extraction
- Partial Validation


Why LCMSMS?

- Why an LCMSMS based assay?
 - ELISA assays not practical for discovery, no antibodies available yet
 - Challenges with ELISA assays
 - time consuming, expensive to develop
 - require separate assay for each peptide
 - limited linear dynamic range
 - Possible cross reactivity
- Benefits of LC/MS/MS for peptides
 - LCMSMS provides single assay for multiple amyloid peptides
 - Broad linear dynamic range
 - Accurate, precise
 - Universal
 - Faster, cheaper method development




Specific Challenges in Developing an LCMSMS Assay for Amyloid β Peptides

- Extremely hydrophobic
- Very poor peptide solubility
- High level of aggregation
- High level of non-specific binding (NSB)
- Low MS sensitivity
 - Very large and hydrophobic
 - Form many low abundance fragments
- Specificity in matrix




600-01695338 [RF] © www.visualphotos.com

ESI- MSMS Spectra for Amyloid β 1-42



ESI+ MSMS Spectra for Amyloid β 1-42



Amyloid β Peptide 1-42: Identification of Fragments with BioLynx In Silico Fragmentation Model

Proteins - B Amyloid 42.pep - [Sequence Ions]

File Edit Display Search Primary Secondary Window Help

a b c⁺ x y⁺ z

Fragment choice is 4+ b 40 ion


B Amyloid 42
Average Mass = 4513.1183, Monoisotopic Mass = 4510.2855
Residues: 1-42
N-Terminus = H, C-Terminus = NH₂
Fragment ions: Monoisotopic/Average (1000) m/z ratios with 4 positive charge(s).

a	722.8	744.6	773.1	805.1	819.4	837.2	865.4	893.7	908.0	936.3	969.1	993.9	1008.1	1022.4	1047.2	1072.0	1100.3	-
b	729.8	751.6	780.1	812.1	826.4	844.2	872.5	900.7	915.0	943.3	976.1	1000.9	1015.1	1029.4	1054.2	1079.0	1107.3	-
c ⁺	734.0	755.8	784.3	816.4	830.7	848.4	876.7	905.0	919.3	947.6	980.4	1005.1	1019.4	1033.7	1058.4	1083.2	1111.5	-
i	30.0	60.0	87.1	101.1	30.0	44.1	86.1	86.1	30.0	86.1	104.1	72.1	30.0	30.0	72.1	72.1	86.1	44.1
	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42
	Gly	Ser	Asn	Lys	Gly	Ala	Ile	Ile	Gly	Leu	Met	Val	Gly	Gly	Val	Val	Ile	Ala
	18	17	16	15	14	13	12	11	10	9	8	7	6	5	4	3	2	1
x	421.3	407.0	385.2	356.7	324.7	310.4	292.6	264.3	235.9	221.6	193.4	160.6	135.8	121.6	107.3	82.6	57.8	29.5
y ⁺	414.8	400.5	378.7	350.2	318.2	303.9	286.1	257.8	229.4	215.1	186.9	154.1	129.3	115.1	100.8	76.1	51.3	23.0
z	410.5	396.2	374.5	346.0	313.9	299.6	281.9	253.6	225.1	210.9	182.6	149.9	125.1	110.8	96.6	71.8	47.0	18.8


Chain 1 of 1

Positive Ion versus Negative Ion Detection: Specificity 1-42 in Human CSF


250pgmL_hCSF_032910_002x




250pgmL_hCSF_032910_002x



Positive Ion versus Negative Ion Detection: Specificity 1-38 in Human Plasma



Effect of Using Higher m/z Precursor/Fragment: plasma



Xevo™ TQ-S MS operated in ESI+ mode MRM Transitions:

Peptide Name	Precursor Ion 4+	Product Ion 4+	Product Ion ID	Cone voltage (V)	Collision energy (eV)
Amyloid β 1-38	1033.5	1000.3	b 36	33	23
Amyloid β 1-38 N15 IS	1046	1012.5		30	22
Amyloid β 1-40	1083	1053.6	b 39	33	25
Amyloid β 1-40 N15 IS	1096	1066.5		35	22
Amyloid β 1-42	1129	1078.5	b 40	28	30
Amyloid β 1-42 N15 IS	1142.5	1091.5		35	28



Solvent Standard Comparison:
Area counts for 500 pg/mL standard


	Xevo TQ	Xevo TQ-S
amyloid β 1-38	108	1468
amyloid β 1-40	133	1741
amyloid β 1-42	54	1589

Xevo TQ-S

LOD: 12.5 pg/mL solvent standard


in 50/50 elution solvent/water

12pt5pgmL_abetas_072610_001 Sm (SG, 3x5)




Amyloid β 1-42
Peak Area 125

12pt5pgmL_abetas_072610_001 Sm (SG, 3x5)



Amyloid β 1-40
Peak Area 98

12pt5pgmL_abetas_072610_001 Sm (SG, 3x5)



Amyloid β 1-38
Peak Area 130

10 μ L injection

- Introduction
- Mass Spectrometry
- ■ Liquid Chromatography
- Choice of Standard Curve Matrix
- Solid Phase Extraction
- Partial Validation

UPLC® Conditions: Final Method

Column: 2.1 X 150mm, ACQUITY BEH C18 300Å, 1.7µm
Mobile phase A: 0.3% NH₄OH by volume, or 0.1% absolute
Mobile Phase B: 90/10 ACN/mobile phase A
Temperature: 50° C
SNW: 65/25/10 ACN/water/NH₄OH
SNW volume: 600 µL
WNW: 90/10 water/ACN + 0.3% NH₄OH
WNW volume: 600 µL
Flow rate: 0.2 mL/min
Injection mode: partial loop
Injection Volume: 10 µL
Injection Solvent for standards:
SPE elution solvent diluted with water;
SPE elution solvent =
75% ACN, 15% water, 10% NH₄OH


Gradient Table

The screenshot shows the 'Modify ACQUITY Binary Solvent Manager Instrument Method' dialog box. The 'General' tab is active. The 'Solvents' section shows Solvent A2 as Ammonium hydroxide and Solvent B2 as *Acetonitrile. The 'Pressure Limits' section shows Low: 0 psi and High: 15000 psi. The 'Seal Wash' is set to 5.0 min. The 'Run Time' is 8.50 min. The 'Gradient' table is as follows:

	Time (min)	Flow (mL/min)	%A	%B	Curve
1	Initial	0.200	90.0	10.0	Initial
2	1.00	0.200	90.0	10.0	6
3	6.50	0.200	55.0	45.0	6
4	6.70	0.200	55.0	45.0	6
5	7.00	0.200	90.0	10.0	6
6					

At the bottom, there is a 'Comment:' field and 'OK' and 'Cancel' buttons.

Representative Chromatogram: Extracted Spiked Artificial CSF




Representative Chromatography: Extracted Human CSF Sample

Basal Levels in Human CSF


hCSF_pool1_basal_040210_001

1-42 N15 IS




hCSF_pool1_basal_040210_001

1-42




hCSF_pool1_basal_040210_001

1-40 N15 IS




hCSF_pool1_basal_040210_001

1-40




hCSF_pool1_basal_040210_001

1-38 N15 IS



hCSF_pool1_basal_040210_001

1-38



- Introduction
- Mass Spectrometry
- Liquid Chromatography
- ■ Choice of Standard Curve Matrix
- Solid Phase Extraction
- Partial Validation


Average Basal Levels in Human CSF: Comparison of Results with Different Standard Curve Preparations

Standard curves prepared from either artificial CSF or by over-spiking human CSF

Amyloid β 1-38	Conc.(ng/mL) Source 1	Conc. (ng/mL) Source 2
Conc from curve in art CSF	0.82	0.826
Conc from curve in human CSF	0.679	n/a
average	0.750	
std deviation	0.100	
RSD	13.302	

Amyloid β 1-40	Conc.(ng/mL) Source 1	Conc. (ng/mL) Source 2
Conc from curve in art CSF	3.713	3.136
Conc from curve in human CSF	4.335	n/a
average	4.024	
std deviation	0.440	
RSD	10.930	

Amyloid β 1-42	Conc.(ng/mL) Source 1	Conc. (ng/mL) Source 2
Conc from curve in art CSF	1.105	0.763
Conc from curve in human CSF	1.189	0.804
average	1.147	0.784
std deviation	0.059	0.029
RSD	5.178	3.700



*Human CSF corrected for basal level

- Introduction
- Mass Spectrometry
- Liquid Chromatography
- Choice of Standard Curve Matrix
- ■ Solid Phase Extraction
- Partial Validation

Sample Pretreatment

Spike human or artificial CSF + 5% rat plasma samples, mix

Equilibrate samples at room temperature for 30 minutes

Remove 50 μL spiked artificial CSF (to which 5% rat plasma was added) for standard curves, or human/monkey CSF (basal level or over-spike QC's)

Add 50 μL 5M guanidine HCl

Shake samples at room temperature for 45 minutes

Add 50 μL 4% H_3PO_4 in water, mix

Note: 4X less sample is used than TQ MS method: 50 μL instead of 200 μL

Sample Extraction: Final Method

Oasis® MCX μ Elution Plate

Condition: 200 μ L MeOH

Equilibrate: 200 μ L 4% H_3PO_4

Load: 150 μ L diluted sample (pretreated sample: 50 μ L human CSF, 50 μ L 5M guanidine HCl, 50 μ L 4% H_3PO_4 in water)

Wash 1: 200 μ L 4% H_3PO_4

Wash 2: 200 μ L 10% ACN in water (by volume)

Elute: 2 X 25 μ L 75/10/15 ACN/conc. NH_4OH /water (by volume)

Dilute: 25 μ L water

Inject: 10 μ L

Amyloid β Peptide	% SPE Recovery
1-38	94%
1-40	92%
1-42	90%



Key Attributes for High Extraction Recovery

- Wash with no more than 10% ACN
 - Removes polar interferences
 - Does not impact recovery of 1-38 (earliest eluting/least hydrophobic of these 3 amyloid β peptides)
- Elute with no less than 75% ACN
 - Provides required elutropic strength for 1-40 and 1-42 (the more hydrophobic of these 3 peptides)
- Elute with no less than 10% NH_4OH
 - Provides required *solubility* to fully elute 1-42 (least soluble/most hydrophobic of these 3 peptides)
- Addition of rat plasma to artificial CSF to eliminate non-specific binding increases SPE recovery in this matrix from 60% on average to >90%
- Guanidine HCl denaturation
 - Improves reproducibility of method; eliminates aggregation and protein binding

Impact of MS Platform: Standard Curve and QC Range Comparison

200 μ L sample

50 μ L sample

	Xevo TQ	Xevo TQ-S
Standard Curve Range	0.1 to 10 ng/mL	0.025 or 0.05 to 10 ng/mL
QC Range	0.2 to 6 ng/mL	0.04 to 6 ng/mL

- Introduction
- Mass Spectrometry
- Liquid Chromatography
- Choice of Standard Curve Matrix
- Solid Phase Extraction
- ■ Partial Validation


- Standard curves prepared in artificial CSF + 5% Rat plasma
 - Calibration points are 0.025, 0.05, 0.1, 0.25, 0.35, 0.5, 0.75, 1, 5, 7.5, and 10 ng/mL
- Basal levels in human CSF
 - 3 individual sources pooled human CSF samples, 1 source cynomolgous monkey
 - 6 replicates from each source
- QC samples
 - Prepared in each of the human CSF matrices
 - 7 concentrations: 0.04, 0.075, 0.15, 0.2, 0.8, 2 and 6 ng/mL over-spike
 - N=3 for each concentration, from each source of CSF

Representative Chromatogram: Basal Levels of Amyloid β 1-38 in Human and Monkey CSF

Basal Level Amyloid β 1-38


hCSF_pool1_basal_040210_001

MRM of 6 Channels ES+
1033.5 > 1000.3 (a beta 1-38)
2.31e4




hCSF_pool2_basal_040210_001

MRM of 6 Channels ES+
1033.5 > 1000.3 (a beta 1-38)
2.96e4




hCSF_pool3_basal_040210_001

MRM of 6 Channels ES+
1033.5 > 1000.3 (a beta 1-38)
1.33e4



cyno_CSF_pool_basal_040210_001


MRM of 6 Channels ES+
1033.5 > 1000.3 (a beta 1-38)
2.43e4




Representative Chromatogram: Basal Levels of Amyloid β 1-40 in Human and Monkey CSF

Basal Level Amyloid β 1-40


hCSF_pool1_basal_040210_001




hCSF_pool2_basal_040210_001



hCSF_pool3_basal_040210_001




cyno_CSF_pool_basal_040210_001



Representative Chromatogram: Basal Levels of Amyloid β 1-42 in Human and Monkey CSF

Basal Level Amyloid β 1-42



Baseline Levels of Amyloid β in Pooled Human and Monkey CSF

Amyloid Beta 1-38

Replicate #	Human CSF	Human CSF Pool 2 ng/mL	Human CSF	Cyno CSF Pool 1 ng/mL
	Pool 1 ng/mL		Pool 3 ng/mL	
1	1.585	2.354	1.014	1.713
2	1.650	2.103	1.371	1.605
3	1.614	2.464	0.950	1.947
4	1.657	1.939	1.608	1.541
5	1.820	2.158	1.471	1.675
6	1.486	1.995	1.167	1.644
Mean	1.635	2.169	1.264	1.688
Std. Deviation	0.110	0.204	0.262	0.140
% CV	6.7	9.4	20.7	8.3

Amyloid Beta 1-42

Replicate #	Human CSF	Human CSF Pool 2 ng/mL	Human CSF	Cyno CSF Pool 1 ng/mL
	Pool 1 ng/mL		Pool 3 ng/mL	
1	0.519	0.616	0.421	0.675
2	0.421	0.656	0.481	0.621
3	0.542	0.644	0.534	0.623
4	0.471	0.567	0.348	0.659
5	0.476	0.573	0.487	0.700
6	0.561	0.713	0.510	0.688
Mean	0.498	0.628	0.463	0.661
Std. Deviation	0.052	0.055	0.068	0.033
% CV	10.4	8.7	14.7	5.1

Amyloid Beta 1-40

Replicate #	Human CSF	Human CSF Pool 2 ng/mL	Human CSF	Cyno CSF Pool 1 ng/mL
	Pool 1 ng/mL		Pool 3 ng/mL	
1	3.083	4.031	2.541	3.699
2	3.391	3.776	2.593	3.989
3	3.292	3.598	2.580	3.525
4	2.884	3.533	2.612	3.956
5	3.131	3.230	2.508	3.284
6	3.656	3.619	2.490	3.595
Mean	3.240	3.631	2.554	3.675
Std. Deviation	0.269	0.266	0.049	0.268
% CV	8.3	7.3	1.9	7.3

Representative Standard Curve: Amyloid β 1-42



Compound name: Amyloid Beta 1-42

Correlation coefficient: $r = 0.998799$, $r^2 = 0.997600$

Calibration curve: $0.887923 * x + 0.0130161$

Response type: Internal Std (Ref 2), Area * (IS Conc. / IS Area)

Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



Representative Standard Curve: Amyloid β 1-42

Name	Type	Std. Conc	RT	Area	IS Area	Response	Conc.	%Dev
blank artificial CSF			5.73	19.7	7.0			
50 pg/mL artificial CSF	Standard	0.05	5.71	230.4	3620.5	0.064	0.057	14
100 pg/mL artificial CSF	Standard	0.1	5.71	390.8	3585.1	0.109	0.108	8.1
250 pg/mL artificial CSF	Standard	0.25	5.71	778.3	3737.3	0.208	0.220	-12
350 pg/mL artificial CSF	Standard	0.35	5.71	1267.3	3693.8	0.343	0.372	6.2
500 pg/mL artificial CSF	Standard	0.5	5.71	1494.7	3566.8	0.419	0.457	-8.5
750 pg/mL artificial CSF	Standard	0.75	5.71	2733.5	4152.0	0.658	0.727	-3.1
1 ng/mL artificial CSF	Standard	1	5.71	3166.8	3792.5	0.835	0.926	-7.4
5 ng/mL artificial CSF	Standard	5	5.72	14773.9	3148.3	4.693	5.270	5.4
7.5 ng/mL artificial CSF	Standard	7.5	5.72	24576.9	3877.0	6.339	7.125	-5
10 ng/mL artificial CSF	Standard	10	5.72	33343.3	3662.5	9.104	10.238	2.4

Average Deviation Values for all Overspike QC Samples

	QC 0.04 ng/mL	QC 0.075 ng/mL	QC 0.15 ng/mL	QC 0.2 ng/mL	QC 0.8ng/mL	QC 2 ng/mL	QC 6 ng/mL
Amyloid β 1-38 Human CSF 1 and 2	2.3	5.8	-3.2	7.3	14.8	5.1	13.1
Amyloid β 1-40 Human CSF 1 and 2	-0.8	-3.2	-1.9	2.5	-2.6	-4.2	-3.8
Amyloid β 1-42 Human CSF 1 and 2	1.3	13.4	-3.6	5.6	2.0	-0.6	-0.2

- **Single flexible LC/MS/MS platform developed for simultaneous quantification of multiple amyloid peptides**
- Highly selective sample preparation based on mixed-mode SPE
- Improved MS selectivity using positive ion mode and sequence ion fragments
- High sensitivity using new MS platform
- Highly reproducible, accurate, and precise
- Sample pretreatment and choice of SPE and LC solutions eliminate handling (NSB, losses, etc.) problems
- One method for multiple A β peptides, advantage over ELISA which requires individual assays for each peptide
- Fast, simple sample prep

Total Solution for Regulated Bioanalysis - Bioanalysis of Peptides

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Sample Prep

- Kits and protocols
- μ Elution plate
- No evaporation & reconstitution
- 15x concentration



UNIFI

- IntelliStart and BioLynx functionalities
- Multiple charge precursor ions

ACQUITY UPLC I Class & Columns

- Lower dispersion
- Dedicated PST columns
- Longer columns if needed
- Lowest carryover



Xevo TQ-S

- Stepwave
- Mass range
- Minimal dwell time
- RADAR
- PIC



SERVICE

- Application support and troubleshooting



- Mary E. Lame, Neuroscience Research Unit, Pfizer Global R&D