

# **Validated bioanalysis for therapeutic antibodies by LC-MS: Fab-selective proteolysis nSMOL**


**SHIMADZU Bioscience Research Partnership (SBRP),  
Shimadzu Scientific Instruments, Bothell, WA, USA**

**Takashi Shimada, Ph.D.  
Noriko Iwamoto, Ph.D.**


**E-mail: [tashimada@shimadzu.com](mailto:tashimada@shimadzu.com)**

# Market trend of antibody drugs

Market size (2017)	2025 forecast
US M\$ 84,500	114,600
Immune checkpoint inhibitor US M\$ 10,566	56,530



# Strategy matching of MS technology to medicine and biologics



# Development of mAb bioanalysis

- **Clinical demands in antibody treatment**
  - ✓ Indicator of drug efficacy
  - ✓ Decision of dosing level
  - ✓ Drug distribution in plasma and tissue
- **Novel bioanalysis for clinical pharmacokinetics**
  - ✓ Structure similarity and sequence specificity for antibody CDR-targeting strategy
  - ✓ Independent of a variety of antibodies
  - ✓ Structure-indicated MS analysis
  - ✓ Clinical PK and discovery for antibody drugs
  - ✓ Regulated LCMS bioanalysis


# FDA Guidance Finalized on May 24, 2018

- Bioanalytical Method Validation Guidance for Industry have been finalized by FDA


## Summary

	Protein LBA	Small molecule LCMS	Protein LCMS
	Same as previously	Almost parameters are same as small molecules.	
<b>Calibration curve</b>		Blank, zero, and six calibration, including LLOQ in every run <b>±15%</b> of theoretical concentration, <b>±20%</b> of LLOQ <b>75%</b> and a minimum of six non-zero calibrator levels should meet the above criteria in each validation run.	
<b>Quality Controls</b>		<b>±15%</b> of theoretical concentrations, <b>±20%</b> at LLOQ for accuracy; and within <b>15%</b> CV, within <b>20%</b> CV for LLOQ	

# Immunoglobulin structure




# Frequency of amino acid substitution in Fv region





# Concept representation of nSMOL

## nano-surface and molecular-orientation limited proteolysis


Trypsin immobilization  
on the surface



Minimizing sample complexity  
into LCMS analysis



Nishio K. Coll Surf B, 2008



CDR Accessible surface

Fab orientation to the solution

IgG collection

Small diameter pore  
100 nm


HOT articles in Analyst

By Katherine Dunn, Publishing Assistant.

Iwamoto N. et.al. Analyst, 2014



# Fv-selective detection by nSMOL




Peptide configuration in Nivolumab 3D structure  
Detected signature peptides from H-chain and L-chain

# Benefit of nSMOL bioanalysis


	nSMOL	Affinity capture	ELISA	
<b>Method R&amp;D</b>	Collection or detection Ab	<b>No</b>	Individual 6-10 months	
	Cross reactivity	<b>No</b>	Yes	
	Effect of ADAs	<b>No</b>	Yes	
	Pre validation	<b>1-3 days</b>	<b>1-3 days</b>	2-3 weeks
<b>Validation</b>	Full validation	3-4 weeks	3-4 weeks	
	Sample prep	<b>Dilution</b> 3-5 hours	Denature, Chemical modification 3-5 hours	<b>Dilution</b> 2-4 hours
<b>Feature</b>	Internal standard	<b>Universal</b>	Individual	
	Dynamic range	Wide	Wide	Narrow
	Selectivity	High	High	Middle-High
	Multiplex assay	<b>Yes</b>	Additional collection Ab	Additional collection /detection Ab

# nSMOL advantage in LCMS bioanalysis


## BV peptide recovery



**BV proteolysis yield: 97%**  
(by SDS-PAGE, densitometry test)



## Calibration curve in BV bioanalysis




# Progress of nSMOL project and clinical trials

Method setting	Full validation	Paper in progress	Clinical trials
25	20	Accepted: 12 In submission: 2	In submission: 3 In trials: 11

Antibodies	Antibodies
Trastuzumab, T-DM1 (0.06-250 µg/ml, Anal Methods, J Pharm Biomed Anal)	Infliximab (0.29-300 µg/ml, Curr Pharm Biotechnol)
Bevacizumab (0.15-300 µg/ml, Drug Metab Pharmacokinet)	Biosimilar of Infliximab (0.29-300 µg/ml, Curr Pharm Biotechnol)
Cetuximab (0.58-300 µg/ml, Bioanalysis)	Etanercept (0.20-100 µg/ml, Pharmacol Res Perspect)
Nivolumab (0.15-250 µg/ml, J Chromatogr B)	Abatacept (0.40-100 µg/ml, Pharmacol Res Perspect)
Rituximab (0.58-300 µg/ml, Biol Pharm Bull)	Tocilizumab (0.78-200 µg/ml, J Pharm Biomed Anal)
Brentuximab vedotin, multiplex assay (0.58-300 µg/ml, Clin Pharma Biopharma)	Coexistence with anti-drug antibodies (Anal Biochem)

# Overall antibody PK for efficacy biomarker

Drug level in disease tissue



Comparative PK  
Dosing strategy

Blood monitoring

Effect of anti-drug antibody  
Cytokine quantitation  
Immune cell profiling

Days

# Acknowledgment

## National Cancer Center

- Akinobu Hamada, Ph.D.
- Hitoshi Nakagama, M.D., Ph.D.



## Kyoto University Hospital

- Kazuo Matsubara, Ph.D.
- Atsushi Yonezawa, Ph.D.



## CMIC Pharma Science Co., Ltd.

- Naoe Yamane, Ph.D.



## Tamagawa Seiki Co., Ltd.

- Ms. Misato Hara
- Mr. Naohiro Hanyu



## Earle A. Chiles Research Institute, Providence Cancer Center

- William Redmond, Ph.D.
- Eric Tran, Ph.D.
- Bernard Fox, Ph.D.
- Hong Ming Hu, Ph.D.
- Yoshinobu Koguchi, M.D., Ph.D

