



**ANAQUANT**

SOLUTIONS FOR BIOANALYSIS

LC-HRMS approach to identify and  
quantitate potential markers for pre-clinical  
and clinical applications in heterogeneous  
samples

Quentin ENJALBERT, PhD

# Context – dermo pharma



- Define potential protein markers for pre-clinical and clinical trials to evaluate the effect of skin treatment
- D-Squame<sup>®</sup> tape samples from different patients (uppermost layer of the skin)
- Main issue: huge heterogeneity in samples
- 6 samples
  - ▀ Control patients
  - ▀ Treated patients



# Objectives & strategies



## Objectives:

- ▀ Define potential protein markers for pre-clinical and clinical trials to evaluate the effect of skin treatment

## Strategies:

- ▀ Use LC-MS for potential biomarker identification
- ▀ Perform a consistent sample comparison that might be applied to higher sample number.



# What can be done in LC-MS?



- DDA or DIA for ID and relative quantification between samples

No standard - No quantitative evaluation – Difficult to compare samples

- Targeted for absolute quantitation

Need ID step first


Best of both worlds?

- Perform ID and quantitative evaluation in the same run



# Top3 strategy in DDA

2 DIFFERENT PROTEINS WITH THE SAME CONCENTRATION IN MATRIX: SUM OF THE 3 BEST FLYER PEPTIDES IS CONSTANT (+/- 10% (SILVA ET AL., 2006))



Universal relative quantification within a sample


Need of an Internal Standard for quantification evaluation



# How to make a universal standard?

🔗 Best universal standard?

- Several proteins over a broad dynamic range




- Complicated to produce 18 calibrated proteins with different concentrations
- Batch-to-batch reproducibility
- Stability unknown

Move to peptide level?



# How to make a universal standard?

- 🔗 Move to the peptide level
  - 🔴 Easier to produce
  - 🔴 3 peptides for each protein: 18 proteins → 54 peptides



How to deal with 54 calibrated peptides over a dynamic range?  
USE OF READYBEADS® TECHNOLOGY

# Workflow with the Universal Quantification READYBEADS



Standards



DIGESTION




DDA


ANALYSIS

Accession Number	Molecular Weight	Protein Grouping Ambiguity
1	28100	64
2	28100	51
3	28100	51
4	28100	51
5	28100	51
6	28100	51
7	28100	51
8	28100	51
9	28100	51
10	28100	51
11	28100	51
12	28100	51
13	28100	51
14	28100	51
15	28100	51
16	28100	51
17	28100	51

Protein IDs



QUANTIFICATION EVALUATION



Identification and quantification evaluation in a single analysis  
Based on classic shotgun workflow



# Identification results



- Sample prep done on 1.5 µg total protein amount
- Injection on nano-LC-Qexactive HF system (Thermo). Column acclaim pepmap, 0.75µm \*50cm.
- 1 hour gradient.

Number of identified proteins in 6 D-Squame® samples						
	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	Sample 6
Number of identified proteins with 2 peptide	161	133	324	365	267	71

Large heterogeneity between samples, even in a same group.  
Difficult to perform classic normalization.





# Quantification evaluation



Independent and individual protein quantification evaluation for each sample

Proteins ID (Confidential)	Sample 1 (µg/mL)	Sample 2 (µg/mL)	Sample 3 (µg/mL)	Sample 4 (µg/mL)	Sample 5 (µg/mL)	Sample 6 (µg/mL)
A	0,002	ND	0,003	0,005	0,002	ND
B	0,008	0,001	0,015	ND	0,002	0,001
C	0,008	ND	0,006	0,005	0,001	ND
D	0,150	0,002	0,029	0,019	0,002	0,002
E	0,010	0,004	0,032	0,019	0,014	ND
F	0,025	0,008	0,045	0,044	0,007	ND
G	ND	0,001	0,016	0,025	0,003	ND
H	0,062	0,009	0,042	0,041	0,002	ND

Direct comparison possible thanks to the quantification evaluation assessment



**ANAQUANT**

SOLUTIONS FOR BIOANALYSIS

# Take-home messages

1. Protein ID and Quantification evaluation within the same LC-MS run.
  - In this case up to 361 proteins ID and with an estimated concentration within the sample
2. Intra and Inter sample quantification evaluation
3. Universal: can be applied to any taxonomy, compatible with DDA and DIA
4. Easy to use workflow
5. Sample comparison time independent

