




Faster Generation of Anti-Drug Antibodies Using SpyTag Technology

Dr. Sarah-Jane Kellmann, Scientist R&D, Bio-Rad Laboratories, Inc

EBF Barcelona, November 25, 2021

In Vitro Antibody Generation

- Generation of monoclonal antibodies for research and in vitro diagnostic use
- HuCAL PLATINUM® phage display technology
- Well-established automated procedures
- Fab antibodies ready in 8 weeks
- Over 17 years of experience
- Generation of almost 50,000 different antibodies




SpyTag

Modular Antibody Assembly Platform


SpyTag Technology

- Based on *Streptococcus pyogenes* fibronectin-binding protein FbaB
- Contains spontaneous isopeptide bond between Lys and Asp




Zakeri B et al. (2012). Peptide tag forming a rapid covalent bond to a protein, through engineering a bacterial adhesin. PNAS 109, E:690-697.

SpyTag-SpyCatcher — Fast, Spontaneous Reaction




- Spontaneous (autocatalytic) reaction
- Covalent isopeptide bond formation, irreversible
- Fast, quantitative reaction
- pH 5–8, temperature +4–37°C
- Robust to buffer conditions, $\text{Ca}^{2+}/\text{Mg}^{2+}$ not needed
- Robust to detergents
- Reaction occurs also inside cells (in vivo)

One Antibody, Multiple Formats in Less Than One Hour




Site-Specific Conjugation

- Cysteines introduced into SpyCatcher
- Site-specific labeling
- Controlled degree of labeling
- High batch-to-batch consistency
- No modification of antibody binding site




SpyTag Fab and SpyCatcher Coupling Reaction




Coupling at room temperature
Ratio Fab : Catcher = 2.5:1

AnyKD Criterion Gel
Nonreduced samples
3 μ g protein loaded per lane



Generation of Anti-Idiotypic Antibodies

Anti-Idiotypic Antibodies



Type 1

- Anti-idiotypic antibody
- Paratope-specific
- Inhibitory
- Neutralizing
- Detects free drug

Type 2

- Anti-idiotypic antibody
- Not paratope-specific
- Not inhibitory
- Detects total drug (free, partially bound, fully bound)

Type 3

- Drug target complex-specific
- Not inhibitory
- Detects bound drug exclusively

A Proven Workflow for Anti-Idiotypic Antibody Generation

Selection strategies to generate binding types 1, 2, and 3

Additional characterization

- Affinity ranking and measurement
- Selection of best pairs

Extended quality control to meet critical reagent requirements

Expertise demonstrated with a portfolio of ready-made anti-biotherapeutic antibodies

Autoimmune disease

- Adalimumab (Humira)
- Certolizumab pegol (Cimzia)
- Etanercept (Enbrel)
- Golimumab (Simponi)
- Infliximab (Remicade)
- Rituximab (Rituxan)
- Tocilizumab (Actemra)
- Vedolizumab (Entyvio)
- Ustekinumab (Stelara)

Age-related macular degeneration

- Ranibizumab (Lucentis)

Allergic asthma

- Omalizumab (Xolair)

Respiratory disease

- Palivizumab (Synagis)

Cancer

- Alemtuzumab (Campath)
- Bevacizumab (Avastin)
- Cetuximab (Erbix)
- Ipilimumab (Yervoy)
- Nivolumab (Opdivo)
- Panitumumab (Vectibix)
- Pembrolizumab (Keytruda)
- Trastuzumab (Herceptin)

Multiple sclerosis

- Natalizumab (Tysabri)

Osteoporosis

- Denosumab (Prolia)

PNH, aHUS (rare diseases)

- Eculizumab (Soliris)

[bio-rad-antibodies.com/biotherapeutics](https://www.bio-rad-antibodies.com/biotherapeutics)

Generation of Anti-Drug Antibodies

- Antibodies specific for biotherapeutic drugs play an important role in preclinical and clinical development
- Used for pharmacokinetic (PK) studies and for the development of immunogenicity assays
- Advantages of integrating SpyTag technology
 - Rapid site-specific labeling
 - Conversion into synthetic IgG in one hour
- Example: Daratumumab (Darzalex)
 - Fully human IgG1 kappa
 - Generated with transgenic mouse technology
 - Developed by Genmab, distributed by Janssen Biotech
 - Treatment of multiple myeloma
 - Binds CD38, causing cells to apoptose

Anti-Daratumumab Antibodies Project — Primary Screening

Compute																									Auto Fill
	Plate:	Plate1	Plate2	Plate3	Plate4	Plate5	Plate6	Plate7	Plate8	Plate9	Plate10														
	Plate ID:	444.12A 10890-1#444.12A 10889-1#444.12A 08060-1#444.12																							


Green: Signal 10-fold above background (73/368, 20%)

Yellow: Signal 5-fold above background


18 clones sequenced → 12 different antibodies

Comparison of Antibody Format Performance in PK Bridging ELISA

Fab-Spy-FcCatcher-HRP





IgG-HRP



- Similar performance of IgG and FcCatcher formats
- Higher sensitivity with FcCatcher due to higher labeling efficiency (site-specific!)

Comparison of Antibody Format Performance in ADA Bridging ELISA



- Performance of anti-daratumumab antibodies compares well across FcCatcher and IgG formats
 - Small scale coupling reaction for many clones possible
 - Rapid conversion of many clones and fast identification of most promising candidates
- Faster, more economical, and more candidates tested than in the conventional workflow

Summary — Recombinant Antibodies with SpyTag Technology

- Fast and robust site-specific labeling of antibodies
- Access to bivalent Fabs or synthetic IgGs within one hour
 - Rapid conversion of many clones, and fast identification of most promising candidates
- Fast and easy switching of species or isotype
- Similar or better performance than corresponding controls in various assays (e.g., ADA, PK)
- Defined product
 - Controlled degree of labeling
 - High batch-to-batch consistency
 - No modification of antibody binding site
- One recombinant Fab gives access to a multitude of experimental setups
 - Choose the best matching Catcher from an “engineering toolbox”



Visit bio-rad-antibodies.com/TrailBlazer for more information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.
HuCAL and HuCAL PLATINUM are trademarks of MorphoSys AG in certain jurisdictions.
© 2021 Bio-Rad Laboratories, Inc.