

Analysis of SARS-CoV-2 using LC-MS Peptide Enrichment for Clinical Research


Dominic Foley
Principal Scientist
Waters, Wilmslow, UK

Overview of SARS-CoV-2 and the role of LC-MS/MS in its detection for clinical research

Overview of a Kit-based LC-MS RUO Workflow for detecting SARS-CoV-2 NCAP Peptides

Performance Characteristics of the Kit-based LC-MS RUO Workflow

SARS-CoV-2 virion and it's proteins





- Spike Protein (S)
 - Recognizes human angiotensin-converting enzyme 2 in the initial stage of infection
- Nucleocapsid Protein (NCAP or N)
 - Involved in (viral) replication and transcription of the genome
 - Structural component of the viral particle
 - Essential for survival

Protein	Intensity	Percentage	MW (kDa)
P0DTC9 NCAP_SARS2	3.1998E+11	88.6	46
P0DTC2 SPIKE_SARS2	17745000000	4.9	141
P0DTC5 VME1_SARS2	13582000000	3.8	25
P0DTD1 R1AB_SARS2	4608000000	1.3	794
P0DTD2 ORF9B_SARS2	2910500000	0.8	11
P0DTC3 AP3A_SARS2	1941200000	0.5	31
P0DTC6 NS6_SARS2	307510000	0.1	7
P0DTC4 VEMP_SARS2	49604000	0.0	8

<https://www.nejm.org/doi/full/10.1056/NEJMci2007042>
<https://www.biorxiv.org/content/10.1101/2020.04.23.057810v1.full.pdf>


Why LC-MS for SARS-CoV-2 Analysis?

- SARS-CoV-2 is a protein-rich viral particle
 - Orthogonal quantitative read-out of digested peptides that correspond to the protein(s) of interest
 - MS based detection compatible with conventional small molecule quantitative systems
- Multi-analyte detection capabilities; other viruses or biomarkers could potentially be detected within the same experiment
- As LC-MS directly detects proteins and peptides, contamination isn't amplified, as might be the case with PCR



```
>sp|P0DTC9|NCAP_SARS2 Nucleoprotein
MSDNGPQNQRNAPRITFGGFSSTGNSQNGERSGARSQRRPQGLPNNNTASWFTALTQHG
KEDLKFFRGQGVPIINTSSPDDQIGYYRATRRIIRGGDGKMDLSPRWYFYLLGTGPEAG
LPYGANKDGIIVVATEGALNTPKDHIGTRNPANNAIVLQLPQGTTLPLKGFYAEGRGGS
QASSRSSRSRNSRNSTPGSSRGTSPPARMAGNGGDAALALLLDRLNQLSKMSGKGGQ
QQGQTVTKSAAEASKKPRQKRTATKYVHVTAQFGRGPEQTQGNFGDQELIRQGTDYKH
WPQIAQFAPSASAFFGMSRIGMEVTPSGTWLTYTGAIKLDKDPNFKDQVILLNRHIDAY
KTFPPTEPKKDKKKKADETOALPQKQKQQTVTLLPAADLDDFSKQLQQSMSSADSTQA
```

Converting Cts to mols for LC-MS




Ct	*Viral Load, RNA/μL	*NCAP Protein amol/μL	*NCAP Peptide amol/μL
20	1.74x10 ⁶	2,900	2,900
25	5.59x10 ⁴	92.8	92.8
30	1.79x10 ³	2.97	2.97
35	57.4	0.0953	0.0953
40	1.84	0.00305	0.00305

*estimates based on cell culture studies

SARS-CoV-2 Workflow Strategies

Enriched and Non-Enriched Workflows





Waters SARS-CoV-2 LC-MS Workflow

Research Use Only

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™


- **Direct detection and quantification** of 3 signature NCAP peptides: ADE, AYN, NPA
- **NCAP proteins:** most abundant and least prone to genetic mutation
- **No replication** of target analytes: reduced contamination and erroneous results



LC-MS/MS Conditions

■ ACQUITY UPLC I-Class FTN

- Column: ACQUITY Premier Peptide BEH C₁₈, 300Å, 1.7µm, 2.1mm x 30mm
- Mobile Phase A: Water with 0.1% formic acid
- Mobile Phase B: Acetonitrile with 0.1% formic acid
- Column Temperature: 40°C
- Injection volume: 20µL
- Run time: 1.8 minutes (2.5 minutes cycle time)



■ Xevo TQ-XS

- Ion mode: ESI+
- Capillary Voltage: 0.5 kV
- Desolvation Temperature: 600°C
- Desolvation Gas: 1000L/Hr
- Source Temperature: 150°C
- Resolution: MS1 & MS2 0.75 FWHM




Peptide	MRM		Cone (V)	Collision (V)
ADE	564.8 > 400.2	Quantifier	35	19
	564.8 > 584.4	Qualifier	35	20
	569.8 > 410.2	SIL	35	19
AYN	563.8 > 679.4	Quantifier	35	19
	563.8 > 892.5	Qualifier	35	19
	568.8 > 689.4	SIL	35	19
NPA	687.4 > 841.5	Quantifier	35	18
	687.4 > 865.5	Qualifier	35	23
	690.4 > 849.5	SIL	35	18
AYE*	588.8 > 400.2	Quantifier	35	19
	588.8 > 777.4	Qualifier	35	15
	569.8 > 410.2	SIL (ADE)	35	19

*Delta variant mutation D377Y

Calibration Performance


- Peptide calibration lines were prepared from 3 – 50,000 amol/μL in VTM (Virus Transport Media)
- Correlation coefficients (r^2) were >0.99

Correlation coefficient: $r = 0.998833$, $r^2 = 0.997667$
Calibration curve: $0.0079962 \cdot x + 0.00127215$
Response type: Internal Std (Ref 4), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Include, Weighting: $1/x^2$, Axis trans: None




ADE

Correlation coefficient: $r = 0.997037$, $r^2 = 0.994082$
Calibration curve: $0.0047293 \cdot x + 0.00113152$
Response type: Internal Std (Ref 5), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Include, Weighting: $1/x^2$, Axis trans: None



AYN


Correlation coefficient: $r = 0.997277$, $r^2 = 0.994560$
Calibration curve: $0.00850259 \cdot x + -0.00154745$
Response type: Internal Std (Ref 6), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Include, Weighting: $1/x^2$, Axis trans: None




NPA

Analytical Sensitivity


- Analytical sensitivity at 3 amol/μL* in VTM after sample enrichment for 3 peptides
 - ADE and AYN provide the greatest analytical sensitivity, with peaks detected in all qualifier traces
 - NPA provides lower analytical sensitivity, with much lower qualifier ion peak areas



ADE



AYN



NPA

*amol= 10⁻¹⁸ mole

Precision

NCAP Protein Derived Peptides and Spiked Peptides



- In-house QCs using spiked NCAP and spiked peptides (ADE, AYN and NPA) into VTM at 3, 10, 400 and 25,000 amol/μL
 - Precision was performed over 5 separate occasions, with 5 replicates at each concentration
 - Spiked Peptides
 - Intra-day and inter-day precision was ≤17.4% CV at 3 amol/μL and ≤12.4% CV from 10 – 25,000 amol/μL
 - Spiked NCAP Protein
 - Intra-day precision was ≤18.8% CV at 3 amol/μL and ≤10.8% CV from 10 – 25,000 amol/μL
 - Inter-day precision was ≤18.8% CV from 3 – 25,000 amol/μL

	Spiked Peptide (amol/μL)							
	Intra-Day Precision				Inter-Day Precision			
	3	10	400	25000	3	10	400	25000
ADE	9.0%	5.0%	2.0%	2.9%	12.9%	8.4%	2.6%	3.3%
AYN	12.5%	6.8%	2.4%	3.0%	15.2%	10.2%	6.8%	4.7%
NPA	17.1%	12.4%	3.0%	4.5%	17.4%	10.9%	5.2%	3.6%

	Spiked NCAP Protein (amol/μL)							
	Intra-Day Precision				Inter-Day Precision			
	3	10	400	25000	3	10	400	25000
ADE	18.4%	7.8%	2.5%	3.1%	11.5%	16.8%	18.8%	11.8%
AYN	13.2%	10.2%	2.4%	2.9%	11.6%	17.6%	18.5%	11.1%
NPA	18.8%	10.8%	4.1%	3.8%	11.6%	15.4%	18.2%	13.3%

Proof of concept


Results based on a similar workflow



Rapid and Sensitive Detection of SARS-CoV-2 Infection Using Quantitative Peptide Enrichment LC-MS/MS Analysis.
Hober et al. <https://elifesciences.org/articles/70843>

The Delta Variant



- As of July 2021, the SARS-CoV-2 Delta (B.1.617.2) variant rapidly became the dominant form of the virus across the globe



- Previous VoCs did not impact the selected peptides ADE, AYN and NPA
- Upon examining the Delta variant, it was noted that the ADE peptide was affected, with increased prevalence of an amino acid substitution adjacent to the N-terminus (D377Y)
 - Original ADE: ADETQALPQR
 - Variant D377Y: A**Y**ETQALPQR

Analysis of Delta Variant NCAP Peptides

- The variant peptide AYE (D377Y) was synthesized and added to a solution containing ADE, AYN and NPA
- All peptides were then added to viral transport medium (VTM) and extracted using the SARS-CoV-2 LC-MS Kit (RUO) and analyzed over a single experimental run
- Skyline was used to identify MRM transitions for the variant peptide and analysis was performed using the ACQUITY UPLC I-Class and Xevo TQ-XS



- No significant changes to the method are required to identify the new peptide target
- The variant peptide can be captured by the SISCAPA ADE wildtype monoclonal antibody
- SARS-CoV-2 LC-MS Kit (RUO) allows the **identification** of Delta variant and its **differentiation** from previous variants of concern in one experiment

A kit-based LC-MS/MS method can be used to directly detect and quantify SARS-CoV-2 NCAP peptides for clinical research

The SARS-CoV-2 LC-MS Kit (RUO) allows direct detection and quantification of three SARS-CoV-2 NCAP peptides with excellent analytical sensitivity, analytical selectivity, reproducibility and speed.

LC-MS can be used as a complementary approach to PCR for research, with the added benefit of using the multi-analyte detection capabilities of LC-MS/MS to also monitor biomarkers as part of SARS-CoV-2 longitudinal research studies

Acknowledgements

- Waters – Chemistry, SciOps, UK Sales
- Cov-MS Consortium
- University Collaborations
- UK NHS
- SISCAPA

SISCAPA®
assay technologies



UNIVERSITY OF
Southampton



viapath

