

Sense and nonsense of a geometric or arithmetic mean medium QC

*Presenter: Peter van Amsterdam,
on behalf of the EBF*

Focus Workshop

(In collaboration with the AAPS and JBF)

**Industry input into ICH M10: Experimental data as the
cornerstone for a science driven bioanalytical guideline**

The Altis Grand Hotel Lisbon,
Portugal September 24-26, 2017

The QC Problem/Challenge/Issue

Quality control (QC) sample (EMA definition)

A spiked sample used to monitor the performance of a bioanalytical method and to assess the integrity and validity of the results of the unknown samples analysed in an individual batch.

Sample Analysis - 3 QC levels:

- QC-LO (within 3xLLOQ; no debate)
- QC-ME (somewhere in some middle; mismatch industry practice & guidance/guideline)
- QC-HI (at least 75% of ULOQ; no debate)

Method Validation – 3 QC levels + LLOQ (+ ULOQ):

- Similar to above (iaw: QC-ME)

What is in the guidelines:

EMA 2011: within three times the LLOQ (low QC), **around 30 - 50%** of the calibration curve range (medium QC), and at least at 75% of the upper calibration curve range (high QC)

ANVISA 2012: amostra de controle de qualidade de média concentração (amostra de CQM): amostra de matriz adicionada do analito em concentração **próxima à média entre os limites inferior e superior** de quantificação
(*Medium-concentration quality control sample (CQM sample): Matrix sample of the analyte in a concentration **close to the mean between the Lower and upper quantification***)

FDA 2013 draft: one within three times the LLOQ (low QC), one in the **midrange** (middle QC), and one approaching the high end (high QC) of the range of the expected study concentrations


MHLW 2012 Chromatography: The low-level should be within 3 times the LLOQ, the **mid-level is around the midpoint** on the calibration curve, and the high-level should be at least 75% of the upper limit of the calibration curve

MHLW 2013 LBA: The low-level should be within 3 times the LLOQ, **the mid-level is near the midpoint** on the calibration curve, and the high-level should be at least one-third of the ULOQ of the calibration curve

Pharmacokinetics: what is important

Conc	C _{max}	AUC _{0-t}	AUC _{0-∞}	t _{1/2}	C _{min}
1					
2					
5					
10					
20					
50					
100					
200					
500					
1000					

Adequately high level of accuracy and precision is important over the entire concentration range !



Just an example of a PK study


Arithm	Geom
0	Blank
	50
	100
	200
	500
	1000
2000	2000
4000	
	5000
6000	
8000	
10000	10000
12000	
14000	
16000	
18000	
20000	20000

qd, bid & tid
SD/MD study,
N=36
2664 samples


Range 50-
20000 pg/ml

Where to
place (how to
distribute) the
Cals?

Arithmetically Distributed Classes




Geometrically Distributed Classes




What we (the EBF community) do:

How do you distribute the concentrations of your calibration standards?

Chromatography




Ligand Binding



- Predominantly use an arithmetical distribution
- Predominantly use a geometrical distribution
- Use both arithmetical and geometrical distributions

Main Activity



- Small molecules/chromatography
- Large molecules/LBA
- Both chromatography and LBA

N=37, 17 European portfolio & 20 Global portfolio

Just an example of a PK study

qd, bid & tid SD/MD study,
N=36, 2664 samples

Range 50-20000 pg/ml

Where to put the QCs?


QC-LO: 100-150 pg/ml

QC-HI: 15000-16000 pg/ml


QC-ME:

- 1000 pg/ml (geomean)
- or 6000 pg/ml (30%)
- or 10000 pg/ml (50% / mean')
- or

Arithmetically Distributed Classes




Geometrically Distributed Classes




Mean, Midpoint, Midrange, 30-50%

Range 1-1000, 10 Cals

Calibration Standards	Mean	Midpoint	Midrange	30-50%
Arithmetically Distributed	500	500	500	300-500
Geometrically Distributed	501	20-50	20-50 or 500	300-500



Arithmetically Distrib Cals			Geometrically Distrib Cals		
Cal	QC-Arithm	QC-Geo	Cal	QC-Arithm	QC-Geo
0			Blank		
100			1		
200	QC-LO	QC-LO	2	QC-LO	QC-LO
300		QC-ME	5		
400			10		
500	QC-ME		20		QC-ME
600			50		
700			100		
800	QC-HI	QC-HI	200		
900			500	QC-ME	QC-HI
1000			1000	QC-HI	




The error distribution

Fitting of calibrator response data to concentration is usually done using the 'Least Squares' approximation. The best fit in the least-squares sense minimizes the sum of squared residuals (the difference between an observed value, and the fitted value provided by a model).


A prerequisite is that the variance in response (y) is similar for each of the concentrations (x). Log transformation e.g. or weighting is thus often necessary for large x ranges.

The parameters of the model, e.g. a and b in $y = a + bx$, are estimates and have a random error and therefore have a confidence interval (CI).



$$s_{\hat{\alpha}} = s_{\hat{\beta}} \sqrt{\frac{1}{n} \sum_{i=1}^n x_i^2} = \sqrt{\frac{1}{n(n-2)} \left(\sum_{j=1}^n \hat{\epsilon}_j^2 \right) \frac{\sum_{i=1}^n x_i^2}{\sum_{i=1}^n (x_i - \bar{x})^2}}$$

$$s_{\hat{\beta}} = \sqrt{\frac{\frac{1}{n-2} \sum_{i=1}^n \hat{\epsilon}_i^2}{\sum_{i=1}^n (x_i - \bar{x})^2}}$$



Where to place your QCs

Confidence interval depends on the distribution of the calibrators. For both arithmetically spaced as well as weighted or log-transformed geometrically spaced CALs the CI is the widest at the extremes.




It makes perfect sense to have a QC-LO at the lower end and a QC-HI at the higher end of the curve.

... the CI is the narrowest at the (transformed or weighted) center of gravity the curve (x-mean,y-mean)....




There is little added value in having a QC-ME in the curve. However smart placing can help detecting lack of fit.


Simulated Data - Arithmetically Distributed Calibration Standards



Simulated Data - Geometrically Distributed Calibration Standards






Simulated Data - Arithmetically Distributed Calibration Standards





What we (the EBF community) do:

Where do you place the QC-ME relative to the calibration standards?



Main Activity





N=37, 17 European portfolio & 20 Global portfolio

Zooming in: what we do

Predominantly a QC-ME at around the geomean of the range of geometrically distributed calibration standards spanning about 3 orders of magnitude.

A QC-ME at around 50% is occasionally used for shorter ranges and/or arithmetically distributed calibration standards.



N=37, 17 European portfolio & 20 Global portfolio



Stating & Concluding that:

- The purpose of fitted calibration curves is to relate response to concentration
- The distribution and range of calibration standards can best be similar to the range and distribution of expected sample concentrations
- The purpose of QCs is to control/manage the quality of analysis batches
- Low and High level QCs express the quality of the assay whereas a Mid level QC could express the suitability of the calibration model provided it is placed wisely along the range

Implies that for the most commonly used geometrically distribute calibration standards:

- A mid level QC at 50% of the highest calibrator gives no extra information on assay performance and will usually always pass and is therefor a waste of effort and materials
- A mid level QC at 50% of the highest calibrator may incidentally give a false positive outcome (4:6:15 / 4:6:20)
- Ambiguity in number and placement of QCs (some labs use 5 levels in production batches) leads to increased costs, timelines and unethical increased use of blank matrices while it is unrelated to the safety and efficacy of medicines

Recommendation

- Geometrically distributed calibration standards: place QC-ME at about the geomean of the lowest and highest standard
- Arithmetically distributed calibration standards: place QC-ME at about the mean of the lowest and highest standard
- Distribute your calibration standards in a geometrically fashion for large(r) concentration spans and arithmetically for short(er) ranges

Arithmetically Distrib Cals			Geometrically Distrib Cals		
Cal	QC-Arithm	QC-Geo	Cal	QC-Arithm	QC-Geo
0			Blank		
10			1		
20	QC-LO	QC-LO	2	QC-LO	QC-LO
30		QC-ME	5		
40			10		
50	QC-ME		20		QC-ME
60			50		
70			100		
80	QC-HI	QC-HI	200		
90			500	QC-ME	QC-HI
100			1000	QC-HI	

Acknowledgement

Abbott, AbbVie, ABL, ABS Labs, Actelion, Almirall, Amanda Wilson, Aptuit, Arcinova, Astellas, AstraZeneca, Bayer, BioAgilytix, Boehringer, Brigitte Buur Rasmussen, Brigitte Buscher, Comac, Covance, CRL, Daniel Mascher, Daniella Stoelner, David Egging, Edwin Hooijschuur, Eginhard Schick, Ferring, Gerhad Paul, Grunenthal, Hans Mulder, Ian Skitt, Janssen, Jo Goodman, Johannes Stanta, Josep M Jansat, Katrin Pöschko, Katrin Schroeter, Klaus Pusecker, LEO, Louis Christodoulou, Luca Ferrari, Lundbeck, Magnus Knutsson, Marco Michi, Marianne Scheel Fjording, Marie-Helene Pascual, Marion Kranenborg, Mark Jean Gnoth, Mellony Dröge, Merck, Mira Doig, Morten Fuch Carlsen, Morten Rohde, MSD, Nathalie Mokrzycki, Novartis, NovImmune, Novo Nordisk, Orion, Pascal Delrat, pharm-analyt, Philip Timmerman, PRA, Prolytic, QPS, Robert Nelson, Roberto Boggio, Roche, Sanofi, Servier, Sirpa Laakso, Stefan Bleich, Steve Williams, Stuart McDougall, Susanne Globig, Synthon, Theo de Boer, Triskelion, Tsvetelina Ivanova, UCB, Ulrich Kunz





Contact: info@europeanbioanalysisforum.eu