

## Development of SARS-CoV-2 ECL Serology Assays

PPD  
Catherine Vrentas, Ph.D.  
Janine Micheli, Ph.D  
Chris Wilson  
Marina Matveyeva


**HELPING DELIVER LIFE**  
**CHANGING THERAPIES**

**PPD**<sup>®</sup>

# SARS-CoV-2 IgG/IgM/IgA Serology Assays

## Assay Format


- + Assays aimed to quantitate the levels of anti-RBD antibodies, specific for SARS-CoV2 infection
- + Utilized commercial RBD reagents and in-house prepared labeled reagents (SulfoTag detection Ab)
- + Unique elements to the approach:
  - Use of MSD/ECL platform
  - Aim to eliminate need for titering
  - Use mix of ADA and biomarker approaches



# Optimization of Assay Format: Need for Custom Normalization

- + Initial assessment revealed that assay potential varied with the subset of COVID-negative serum samples used
- + For example, normal serum samples typically had signals in the thousands, but certain individuals showed very elevated signals (false positive)
- + Optimized buffer composition (base buffer formulation utilizes Blocker Casein), which improved background, but issues persisted

ECL Signal: Negative serum




# Optimization of Assay Format: Need for Custom Normalization



- + The background was not specific to coated wells!
- + Similar background profiles observed for coated wells vs. wells that were blocked only.
- + Therefore, not due to the issue of cross-reactivity from other anti-coronavirus antibodies.
- + Allows for a custom normalization strategy: subtract the signal from blank wells from the coated wells.

## ECL Signal: Negative serum

Black bars: Coated & blocked wells  
Gray bars: Blocked wells




# Advantages of normalization apparent in patient sample comparison



**Red lines = Positive Samples**  
**Blue = Negative Samples**

# Cut point approach demonstrates potential of the method

- + Assessed the performance of the method with background correction across a set of 60 samples, run across 3 days (commercial serum)
- + All samples were COVID-negative: collected in U.S., before Jan. 2020
- + Upper panel: Demonstrates reproducibility of raw (prior to correction) signals
- + Lower panel: Demonstrates utility of the background correction approach.




# Anti-SARS-CoV-2 IgG Serology Assay

## Sample results at 1:20 dilution


- + A panel of negative samples (n = 59) was used to estimate a statistical cut point to discriminate between positive and negative samples
- + This allows an estimate of clinical sensitivity and specificity
- + Samples with PCR confirmed COVID-19 diagnoses (n=19) were used to verify cut point performance

	Negative Sample	Positive Sample
Positive Result	1/59 (1.7%)	18/19 (94.7%)
Negative Result	58/59 (98.3%)	1/19 (5.3%)
Diagnostic Sensitivity	94.7%	
Diagnostic Specificity	98.5%	



# Spiked samples further demonstrate potential of the method

- + Artificially-spiked sample further demonstrate the utility of the method
- + Used an anti-RBD IgG PC spiked into COVID-negative serum samples
- + Even at a low dilution (MRD20), results clustered across numerous individuals
- + Readily able to detect only 200 ng/ml PC antibody vs. negative samples






# Selectivity Validation

- + Positive and negative serum samples were analyzed
- + Titered using a 2-fold dilution scheme
- + 20- to 1280-fold dilutions were tested
- + The calibration curve LLOQ is provided here—can adjust based on desired specificity
- + Samples with results in the quantitation range can be reported as relative concentrations

Statistic	Value	95% CI
Sensitivity	90.0%	55.5% to 99.8%
Specificity	100%	69.2% to 100%
Positive Likelihood Ratio		
Negative Likelihood Ratio	0.10	0.02 to 0.64
Disease prevalence (*)	10.0%	
Positive Predictive Value (*)	100%	
Negative Predictive Value (*)	98.9%	93.3% to 99.8%
Accuracy (*)	99.0%	81.4% to 100%

(\*) These values are dependent on disease prevalence.




**Red lines = Positive Samples**

**Blue = Negative Samples**

# Sample Quantitation: Use biomarker-like approach, reference curve

- + While an ADA cut point approach could be used, observed some differences in signal across plates, but not necessarily tracking with “noise”
- + Therefore, instead of targeting S/N semi-quantitation, used a reference calibration curve on plates
- + Reference curve: screen matrix to remove outliers, spike commercial PC (anti-RBD PC) into matrix
- + Quantitate normalized sample signals from the curve



# Formal validation supports the robustness of the method

- + Inter-assay precision of calibrators (back-calculated values): 2-7%
- + %CV of calibrator raw responses across 17 runs (4 analysts) = 12-18%
- + Run acceptance based on suitability of calibrators (%CV) and adherence to reference ranges as established in validation (similar to ADA HPC, LPC), as well as blank assessment (uses blank matrix pool)
- + No interference from sample lipemia or hemolysis

## Parallelism analysis: increased granularity of quantitation

- + Titer assessment reports within 2-fold range on each side of mean
- + Quantitative approach allows for further granularity of reporting
- + Assessed parallelism of set of samples—dilute into sample diluent
- + Quantitated as diluted into buffer; established against the matrix pool curve (which is spiked into buffer to create the MRD40)


Sample	Dil 20	Dil 40	Dil 80	Dil 160	Dil 320	Dil 640	Dil 1280
1	0	-9.2	< LLOQ	< LLOQ	< LLOQ	< LLOQ	< LLOQ
2	> ULOQ	> ULOQ	0	-3.7	-13.4	-6.4	-12.8
3	0	-29.4	< LLOQ	< LLOQ	< LLOQ	< LLOQ	< LLOQ
4	> ULOQ	> ULOQ	rcv	0	5.7	-6.7	-15.5
5	0	0.30	rcv	9.3	18.0	26.1	< LLOQ
6	> ULOQ	> ULOQ	0	5.3	2.4	-0.96	-6.7
7	> ULOQ	> ULOQ	> ULOQ	0	-9.2	3.5	< LLOQ
8	rcv	0	8.8	27.7	40.3	< LLOQ	< LLOQ
9	> ULOQ	0	-0.85	-7.6	-15.3	-20.8	-19.2



# Assay components have also been adapted to plate-based Nab

## Nab Format:

- Standard cell-based assay examines ability of anti-RBD antibodies in serum to block (pseudo)viral uptake to cells—by preventing ACE2 binding
- Plate-based format assesses ability of serum Nabs to block ACE2-RBD binding interaction
- Used cut point approach to assess presence of neutralizing antibodies in the serum



Copyright, 2020 by Pharmaceutical Product Development, LLC ("PPD"). All rights reserved. This presentation, including the information contained herein and commentary associated herewith ("materials"), is provided as a service of PPD. These Materials, based on publicly available information as well as the knowledge and experience of PPD's employees, have been prepared solely for informational and educational purposes and should not be relied upon for professional advice. Any further use of these Materials requires the express written consent of PPD.