



VNiVERSiDAD
D SALAMANCA



UNIMORE
UNIVERSITÀ DEGLI STUDI DI
MODENA E REGGIO EMILIA

New Analytical Tools for the Study of Bioactive Polyphenols in *Prunus avium* L. (sweet cherry) and Evaluation of their Anti-oxidant Activity on an *In vivo* Model

Dr. Virginia Brighenti, Ph.D.

*Department of Life Sciences, University of Modena and Reggio
Emilia*

6th EBF YSS
24-25th September 2020

Summary

1. Background and aim

2. Analytical characterisation

- Conventional analytical techniques (HPLC-UV-ESI-MS and MS²)
- Electrochemical analysis (sensors)

3. *In vivo* assessment of the anti-oxidant activity on *C. elegans* model

In collaboration with the research group of Prof. Celestino Santos-Buelga of the University of Salamanca

Prunus avium L. I.D.



- **Family:** *Rosaceae*
- **Genus:** *Prunus*
- **Species:** *avium*
- **Nickname:** sweet cherry
- **Fruit:** drupe

Inside sweet cherry fruit




Polyphenols

Sugars

Organic acids

Vitamins

*«Let food be your medicine and medicine be your food»
- Hippocrates*



Gonçalves et al., Sweet Cherry Phenolic Compounds: Identification, Characterization and Health Benefits, in: Atta-ur-Rahman (Ed.), Studies in Natural Products Chemistry, Elsevier Inc., Amsterdam, 2018.

Our goal

- Application of suitable analytical techniques to the **characterization** of sweet cherry polyphenolic fraction **and development of a fast and innovative tool** alternative to conventional techniques.
- Evaluation of the capacity of sweet cherry extract to increase the **resistance to thermal stress** in *C.elegans in vivo* model.



Sample



In this work, cherries belonging to the variety «**Moretta**» were considered. The sample was harvested near the city of Vignola in June 2018 and kept at $-20\text{ }^{\circ}\text{C}$ until the analysis.



Highlights

1. Background and aim

2. Analytical characterisation

- Conventional analytical techniques (HPLC-UV-ESI-MS and MS²)
- Electrochemical analysis (sensors)

3. *In vivo* assessment of the anti-oxidant activity on *C. elegans* model

In collaboration with the research group of Prof. Celestino Santos-Buelga of the University of Salamanca


Workflow

**Extraction of
secondary
metabolites**



**Quantification of the
main classes of
polyphenols
(Uv-Vis)**

**Identification of the
compounds
(HPLC-UV/DAD and
MS and MS²)**

**Electrochemical
analysis**



Spectrophotometric analyses



Which anthocyanins?
Which flavonoids?
Which other phenolics?



Compounds identification

HPLC-ESI-UV/DAD and MS and MS² analysis: combining spectral, MS and fragmentation information.

Chromatographic parameters

Chromatographic apparatus: HPLC Agilent Technologies modular model 1100;

Column tested: Zorbax C₁₈ SB (150 x 4.6 mm, 5 μm, Agilent Technologies),

Mobile phase:

A) H₂O-HCOOH (9:1, v/v),

B) MeOH-H₂O-HCOOH (5:4:1, v/v/v), gradient elution;

Injection volume: 5 μL;

Flow: 1.0 mL/min (for MS and MS² analysis flow was splitted 1:5);


Column temperature: 25°C;




Detection parameters

UV/DAD detection: 280 nm and 320 nm (phenolic acids and flavonoids), 520 nm (anthocyanins);

MS-IT detection: The ion trap mass analyzer was used in the full-scan **positive** (anthocyanins) and **negative** ion modes (flavonoids, phenolic acid) , in the *m/z* range 100-1000. MS² spectra were automatically performed with helium as the collision gas, using the SmartFrag function.




Compounds identification



Peak n°	Anthocyanin
1	Cyanidin-3-O-glucoside
2	Cyanidin-3-O-rutinoside
3	Pelargonidin-3-O-rutinoside
4	Peonidin-3-O-rutinoside

Representative HPLC-UV/DAD chromatogram of **sweet cherry fruit extract** recorded at 520 nm

Compounds identification



Representative HPLC-UV/DAD chromatogram of **sweet cherry fruit extract** recorded at 280 nm

Peak n°	Compound
1	Neochlorogenic acid
2	Catechin
3	3-Coumaroylquinic acid
4	Procyanidin dimer B type
5	5-Caffeoylquinic acid
6	Ferruoylquinic acid
7	cis-4-Coumaroylquinic acid
8	Epicatechin
9	trans-4-Coumaroylquinic acid
10	Procyanidin dimer B type
11	Naringenin-hexoside
12	3,5-Dicaffeoylquinic acid
13	Quercetin-3-O-glucoside
14	Quercetin-3-O-rutinoside
15	3-Coumaroyl-5-caffeoylquinic acid
16	Catechin-hexoside

■ Hydroxycinnamic acids ■ Flavan-3-ols ■ Flavonols ■ Other Flavonoids

Come to the unconventional side!

Conventional techniques

- High costs
- Time consuming sample pre-treatment
- Large sample size
- Large solvents and reagents usage

Electrochemical analysis

- Low costs
- Fast (minimal sample pre-treatment)
- Small sample size
- Limited solvents and reagents usage



**Amperometric
sensors**

Amperometric sensors




Figure. Scheme of a typical circuit used with amperometric sensors.



- **Chemical sensors**
- measure **current response** to detect the **concentration** of an analyte at a fixed potential.
- The **current** is generated by a **RedOx process** between the **sensor surface** and the **analyte**.
- **RedOX reaction** occurs at the **working electrode**, (the real **sensor**) and determines the passage of current which, in particular conditions, is **proportional** to the **analyte concentration**.

**In-line measurement
no sample prep
user-friendly
miniaturization**

Electrochemical conditions


- **Working electrode:** Sonogel Carbon electrode containing Carbon Black (SNGC_CB);
- **Reference electrode:** saturated Ag/AgCl;
- **Auxiliary electrode:** glassy carbon;
- **Working solution:** 0.1 M tartrate buffer solution (TBS), pH= 3.5;

- **Differential Pulse Voltammetry (DPV) parameters:**
 Pulse time: 0.15 s
 Voltage stepsize: 0.6 s
 Modulation amplitude: 0.01 V
 Equilibration time: 5 s




Qualitative analysis of anthocyanins

	compound	E / V
pelargonidin	Pg	0.683
petunidin	Pt	0.441
	Ptg	0.432
peonidin	Pn	0.706
	Png	0.695
delphinidin	Dpr	0.416
	My	0.407
malvidin	Mvg	0.613
	Oe	0.601
cyanidin	Cn	0.473
	Ku	0.487
	Cy	0.508



**Anthocyanins
reference
standard
solution**



**Mashed sweet
cherry fruit**

Summary

1. Background and aim

2. Analytical characterisation of sweet-cherry fruit extracts

3. ***In vivo* assessment of the anti-oxidant activity on *C. elegans* model**

In collaboration with the research group of Prof. Celestino Santos-Buelga of the University of Salamanca

Workflow



**Sweet cherry fruit
(«Moretta» variety)**



**Preparation of the
extract**



***In vivo* assays
on *C. elegans***




**Resistance to
thermal stress**

Preparation of the extract




«Moretta»
sweet cherry
fruit



Dynamic maceration
MeOH/2% HCl
95:5 (v/v)



Concentration



Purification
with water
on C₁₈ column



Final
Cherry Extract
(CE)




Lyophilization




Experimental

- The effects on the **resistance to thermal stress** of the extract obtained from “Moretta” cherries have been evaluated.
- Thermal stress assays were carried out in **wild type N2** of the model nematode *C. elegans* **grown** on NGM agar plates **in the presence and absence** (control group) of the **cherry extract** at different concentrations, on the basis of previous experiments.
- Resistance to thermal stress was evaluated by measuring the **rate of survival**, through the response to a mechanical stimulus, after **6 or 8 hours** of **exposure** to the **lethal temperature of 35 °C** at different stages of the development (2nd and 9th day of adulthood).
- Assays were performed with approximately 100 nematodes for treatment, and each trial was made in triplicate. In each assay, the results were expressed as percentage of living animals.



Results

- **None** of the concentration of the cherry extract tested leaded to **toxic effect** on nematodes.
- The most significant results on the anti-oxidant activity were obtained at a concentration of 25 µg/mL. Indeed, an **increase in survival** in the nematodes treated with cherry extract both on the **second and ninth day** of worm life was evident.
- The most significant result was obtained on the **second day (A)** for a thermal stress of **eight hours** where there was a 25% increase in the survival of treated worms compared to the control. The differences were considered significant at *(p<0,05), **(p<0,01) and ***(p<0,001).



Overall conclusions

- In this work, cherry fruit belonging to the «Moretta» variety were characterized in the composition of their polyphenolic fraction by means of conventional analytical techniques. The **polyphenolic fraction** of «Moretta» cherry fruit encompassed **anthocyanins, hydroxycinnamic acids, flavan-3-ols** and other flavonoids.
- An **innovative analytical tool** based on electrochemical analysis with **SNGC_CB amperometric sensors** for the **fast** determination of anthocyanins in sweet cherry fruit was developed. The sensor proved capable of selectively identify the major class of anthocyanins in cherry fruit, namely cyanidin-3-*O*-glycosides. Future steps will include the quantitative determination of the anthocyanins by using this innovative technique.
- Given the high biological value of the polyphenolic fraction of “Moretta” cherry fruit, its extract was submitted to *in vivo* antioxidant activity assays by using the *C. elegans* model. The extract was able **to significantly increase the adult worm resistance to oxidative and thermal stress** at a concentration of **25 µg/mL**.
- Moretta cherry extract represents a **promising source** of bioactive compounds, thus it might find an application in the nutraceutical industry.

Acknowledgments



Prof. Federica Pellati, Ph.D.
Prof. Davide Bertelli, Ph.D.
Prof. Laura Pigani, Ph.D.
Prof. Chiara Zanardi, Ph. D.
Prof. Stefania Benvenuti, Ph.D.
Irene Laccu
Lucia Marchetti
Corinne Ngnamoko



Prof. Celestino Santos-Buelga, Ph.D.
Prof. Ana Maria González-Paramás, Ph.D.
Francesca Cimurri

FAR 2018 Project: “Sweet cherry and other red fruit as new sources of neuroprotective agents: an interdisciplinary study” funded by UNIMORE.

Thank you!!

...questions?

e-mail contact:

virginia.brighenti@unimore.it

