



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.



Gate that cell!

Requirement for analyzing
Flow Cytometry Data
Reflection on H62
document

Teona Roschupkina



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Outline





- Overview of how data generated
- Poisson distribution for flow
- Variability –accounting for planning
- Gating strategies and their usefulness



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Brief overview of Flow Process

How flow works..







**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Brief overview of Flow Process

	Violet DAPI	Green GFP	Orange Propidium Iodide
Cell 1	200	30	30
Cell 2	230	175	40
Cell 3	250	50	50
Cell 4	180	190	60
Cell 5	200	46	70






**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

**Gating is used to
derive data, but
it helps us clean
up data.**

Gating





**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Variability enemy of statistics

H62 and FDA references to rare cells

“To define a cell population [to gate it] and allow reproducible enumeration it has been recommended that 20-50 clustered events be present...” H62

Desired CV	1%	2.5%	5%	10%	20%
Required Events for the Population of Interest	10 000	1600	400	100	25
Population of Interest Frequency, %	Total Number of Acquired Events				
10	100 000	16 000	4000	1000	250
1	1 000 000	160 000	40 000	10 000	2500
0.1	10 000 000	1 600 000	400 000	100 000	25 000
0.01	100 000 000	16 000 000	4 000 000	1 000 000	250 000
0.001	1 000 000 000	160 000 000	40 000 000	10 000 000	2 500 000



DRUG DEVELOPMENT
SOLUTIONS
Part of Alliance Pharma, Inc.


Poisson distribution
How many cells we
need?!

Poisson counting statistics







$$CV, \% = \frac{\sqrt{N}}{N}$$

$$CV = \frac{\sqrt{11}}{11} = 30.2\%$$







$$CV_{\text{experiment}} = 48.4\%$$



Binomial



Poisson



time



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Variability enemy of statistics

Sources of Variability for cells

- Population separation
- Blood Lysis vs. PBMC
 - Debris and Viability
- Analysts
- Instruments
- Temperature
- Centrifugation, xg force
- Type of assay
 - Intracellular
 - PhosphoFlow
- Logistics
- Sample type
- Analysis of data



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Phosphorylated Biomarker MD

- Very low frequency
- Expressed in cell lines
- Mitosis biomarker
- Predictive in number of tumours as progression marker
- Matrix is rare
- Not every patient expresses biomarker of interest
- This was whole blood assay



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

H62 Specificity and Gating


“...Once gating strategy is determined the results should be verified against reference or comparative method.....”

“...final document should provide how gating was designed to ensure specificity”




**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Spiked healthy individuals



0.177% of total
CV=4.04




0.18% of total
CV=2.71



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.


Spiked healthy individuals

MNC & Spiked Cells



0.03% of total
CV=15.71% of
Blast gate

MNC & Cells





0.03% of Blast gate
CV=8.32%



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Disease Matrix




High Concordance!!!




**DRUG DEVELOPMENT
SOLUTIONS**

Part of Alliance Pharma, Inc.

Disease Matrix



0.1% of blasts
CV=0%



0.01% of blasts
CV=0%

Conclusion



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

- Plan should contain descriptive element of gating: Clustered events, Pattern
- Gating strategy comes first!
- Effective gating strategy for biomarker and disease
- Inclusion of further means to resolve population may improve specificity and sensitivity



**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

Thank you for your attention

Any further questions?

drugdevelopmentsolutions.com



@DDSDrugDev



@drugdevelopmentsolutions



@drugdevelopmentsolutions





**DRUG DEVELOPMENT
SOLUTIONS**
Part of Alliance Pharma, Inc.

**Additional Slide
for question**

**Q: Were there any cases with
Mitotic Marker expressed
aberrantly?**

