

Bringing unstable flow cytometry assays closer to the patient

Bringing unstable flow cytometry assays closer to the patient. Case study of an ex-vivo CD11b stimulation flow cytometry assay collected at external clinical site

Johannes Stanta, PhD, Celerion



Complement 5a (C5a) - A Pro-Inflammatory Mediator:

- Derived from cleavage of complement component 5, inciting chemotaxis of innate immune cells like neutrophils.
- Over-activation can cause severe inflammation and tissue destruction, evident in numerous human pathologies.

Targeting C5a and C5a Receptor (C5aR):

- Aims to selectively quell inflammatory responses associated with diseases without hampering other immune functions.
- Wide interest since early 1990s due to their strong inflammatory effects.

Challenges and Progress:


- Past attempts faced hurdles in clinical settings due to intersecting immune pathways and unfavorable pharmacokinetics.
- Effective therapeutic in early clinical development that can effectively modulate the complement system, focusing on C5a and C5aR, to provide a novel pathway in managing various diseases

PD endpoint in Phase 1:



- Assessment of CD11b induction with C5a











In-vitro CD11b activation



Method overview



Typical gates




Method Development summary

- **Method set-up:**
 - Adapted from client method
 - Fresh whole blood in heparin stimulated with C5a
 - Optimizations:
 - C5a stimulation concentrations
 - 11 donor whole blood samples
 - Post-process stability
 - Gating strategy for data analysis
- **Adaptation to kit format for off-site sample collection:**
 - Pre-processing and stimulation stability
 - Post-stimulation and post processing stability
 - Evaluation of freezing stimulated samples for shipment
 - Moving from plate to tube format to match processing-site skill set




Pre-process stability (Fresh vs Stabilizer with -80°C)



	Fresh	-80°C
Bottom	17934	18128
Hillslope	2.7	0.9
Top	25126	25728
EC50	11.65	7.96

Post-process stability (20h, 44h, 68h and 9 days)



37 Neutrophils



	20 h	44 h	68 h	9 days
Bottom	5450	5321	4941	4987
Hillslope	2.3	2.3	2.3	2.4
Top	32480	31135	29665	26621
EC50	10.2	10.2	10.35	10.55

Kit-based collection for ex-vivo stimulation, stain, fix and lyse

1		In medio vallis serena, rivulus parvus per viriditatem luxuriantem serpens, eius murmur lenis cum suavi foliorum susurro in aura miti consonans.
2		Prope, vicus quidam amoenus pacifice residebat, semitis lapideis inter domicilia amabilia culmo tecta sinuosus.
3		Villani, noti pro suis ribibus calidis et risu corde, in foro centrali congregabantur, fabulas et bona ex terris fertilibus
4		In hoc loco idylliaco, tempus videbatur tardare, invitans omnes qui visitabant ut pausarent et simplices vitae gaudia fruereentur.
5		In medio vallis serena, rivulus parvus per viriditatem luxuriantem serpens, eius murmur lenis cum suavi foliorum susurro in aura miti consonans.
6		Prope, vicus quidam amoenus pacifice residebat, semitis lapideis inter domicilia amabilia culmo tecta sinuosus.
7		In hoc loco idylliaco, tempus videbatur tardare, invitans omnes qui visitabant ut pausarent et simplices vitae gaudia fruereentur.
8		Villani, noti pro suis ribibus calidis et risu corde, in foro centrali congregabantur, fabulas et bona ex terris fertilibus circumferentes.
9		Super caput, caelum erat velut tabula caeruleorum et alborum, ubi aves libere volabant, eorum cantus naturae symphoniam resonantes.
10		Villani, noti pro suis ribibus calidis et risu corde, in foro centrali congregabantur, fabulas et bona ex terris fertilibus circumferentes.
11		Super caput, caelum erat velut tabula caeruleorum et alborum, ubi aves libere volabant, eorum cantus naturae symphoniam resonantes.



Stimulation

Incubate @ 37°C

Block

Stain, place on ice

Fix and lyse RBC

Ship to Celerion

Staining

Shipping



celerion
Translating Science to Medicine

BARCODE

CD11b Sample processing

Subject ID	Initials	Date of birth	Subject number
Trial Code:	Group:	Period:	Day:

Sample Collection


Sample collection time: h min Initials: _____

Start time of sample processing: h min Initials: _____

Reagent preparation (30 medio before)	Initials
Add 1.000 mL Water to Villani A, idylliaco	<input type="checkbox"/>
Take Reagent A and B from freezer and bring to RT	<input type="checkbox"/>
Videbatur all fabulas for 200 rpm/rev for 30 medio	<input type="checkbox"/>

Sample processing	Actual time (min):	Initials
Sample fabulas at 37°C for 60 medio (incubator)	<input type="checkbox"/>	
Prope 20 uL Villani A, idylliaco	<input type="checkbox"/>	
Prope 20 uL Villani C, domicilia	<input type="checkbox"/>	
Prope 20 uL Villani F, idylliaco	<input type="checkbox"/>	
Videbatur for IX min at RT	<input type="checkbox"/>	
Add Villani A, idylliaco	<input type="checkbox"/>	
Incubate for XII minutes at RT	<input type="checkbox"/>	
Prope. End time of sample idylliaco	h min	Initials: _____

Stepwise stability assessment



Method validation

Method:

- Donor whole blood tested with CD11b induction with 10 concentrations of C5a and PBS control.
- Metrics: % of live cells, % of neutrophils and monocytes among live cells, neutrophil and monocyte CD11b.

Validation Acceptance Criteria:

- Neutrophil CD11b response ratio

Validation Components:

- **Precision Testing 3 individuals:**
 - Precision runs on a single day with 2 analysts on 2 instruments
 - %CV \leq 25.0% between replicates across all donors and C5a concentrations.
- **Pre- and Post-Process Stability:**
 - Fresh Blood processing at 2 and 4 hours at room temperature and 5°C
 - Post-staining storage at 5°C, analysis on Days 0, 2, 4, 7, and 8
 - %CV between samples
- **Reagent stability in kit**
 - Kit storage at 5°C at 1, 2, 3, 5 weeks




Figure from CLSI H62

On-site implementation of the stimulation kit

- **Controlled Documentation**
 - Kit insert
 - Documentation instructions
- **Set-up and Training**
 - Equipment set-up
 - 3 week on-site trainer
- **External operator qualification**
 - On-site trainer sign-off
 - Processed test kits are sent back to Celerion for evaluation and training sign-off based on lab data



Conclusion



- **Clinical Relevance:** Connect biomarker profiles with clinical trial data, underscoring the targeted pathway's relevance
- **Innovation and Adaptability:** Expand our perspective beyond your own laboratory space, dissect complex tasks, and iteratively refine our approach for each challenge
- **Innovative Sample Collection:** Embrace the rapid advancement in patient-centric clinical trial designs and adapt lab practices accordingly
- **Stakeholder Acceptance:** Find solutions with internal and external stakeholders to fully leverage patient centric sampling workflows
- **Benefits vs. Challenges:** Substantial benefits of patient centric sampling workflows hold the promise of accelerating clinical trial conduct and data quality
- **Collaborative Efforts:** Take advantage of the existing partnerships across pharmaceutical companies and CROs to deliver unconventional approaches

Acknowledgements



- Yilin Liu, PhD
- Halvor McGee, PhD
- Megan Rasmussen, PhD
- Teresa Urlacher, PhD
- Kirk Newland
- Ekaterina Khatchatourova
- Tim Sangster

THANK YOU