

MORGAN EVANS – SENIOR SCIENTIFIC ADVISOR

Can technology choice make your data SPARCL?

QUALITY

INNOVATION

INTEGRITY

agilexbiolabs.com



AGENDA

01

An introduction to the background literature that inspired this work

02

Tocilizumab Case Studies: Critical Reagents and Methods.

03

Data, Data, Data!

04

Lessons learned and the future perspective.

AGENDA

01

An introduction to the background literature that inspired this work

02

Tocilizumab Case Studies: Critical Reagents and Methods

03

Data, Data, Data!

04

Lessons learned and the future perspective.

AGENDA

01

An introduction to the background literature that inspired this work

02

Tocilizumab Case Studies: Critical Reagents and Methods

03

Data, Data, Data!

04

Lessons learned and the future perspective.

AGENDA

01

An introduction to the background literature that inspired this work

02

Tocilizumab Case Studies: Critical Reagents and Methods

03


Data, Data, Data!

04

Lessons learned and the future perspective.

01

An introduction to the background literature that inspired this work



01 - BACKGROUND

01 - BACKGROUND

Review

For reprint orders, please contact: reprints@future-science.com

Bioanalysis

Evaluation of homogeneous proximity immunoassays for preclinical bioanalysis

Padmanabhan Eangoor*,¹ Sharmistha Das¹ & Vincenzo Pucci¹

¹Pharmacokinetics, Pharmacodynamics & Drug Metabolism, Merck & Co., 33 Avenue Louis Pasteur, Boston, MA 02115, USA

*Author for correspondence: Tel.: +1 617 992 3270; padmanabhan.eangoor@merck.com

Information:

MSD, Gyrolab®, HTRF, AlphaLISA®, and SPARCL™ technologies included.

Evaluation of IgG, IL-8 and Insulin assays.

Parameter	MSD	Gyrolab®	SPARCL™
Assay Time (Hrs)	4-6	1-2	1
Sensitivity	++++	++++	+++
Dynamic Range	3-4	3-4	3-4

01 - BACKGROUND

Review

For reprint orders, please contact: reprints@future-science.com

Bioanalysis

Evaluation of homogeneous proximity immunoassays for preclinical bioanalysis

Padmanabhan Eangoor*,¹ Sharmistha Das¹ & Vincenzo Pucci¹

¹Pharmacokinetics, Pharmacodynamics & Drug Metabolism, Merck & Co., 33 Avenue Louis Pasteur, Boston, MA 02115, USA

*Author for correspondence: Tel.: +1 617 992 3270; padmanabhan.eangoor@merck.com

The opportunity:

Assessment of the MSD, Gyrolab, and SPARCL platforms.


Utilising the same reference material and critical reagents such as capture and detection antibodies.

02



Tocilizumab Case Studies: Critical Reagents and Methods



PK ASSAY FOR THE DETERMINATION OF TOCILIZUMAB IN HUMAN SERUM



PK ASSAY FOR THE DETERMINATION OF TOCILIZUMAB IN HUMAN SERUM



H BioRad Clone, AbD21362





Tocilizumab



BioRad Clone,
AbD21338_hlgG1

PK ASSAY FOR THE DETERMINATION OF TOCILIZUMAB IN HUMAN SERUM



H BioRad Clone, AbD21362

 Tocilizumab

 BioRad Clone, AbD21338_hIgG1





MSD: SulfoTag™

Gyrolab: Biotin & Alexafluor

SPARCL: HRP & Acridan

PK ASSAY FOR THE DETERMINATION OF TOCILIZUMAB IN HUMAN SERUM



H BioRad Clone, AbD21362



Tocilizumab



BioRad Clone, AbD21338_hIgG1



MSD: SulfoTag™

Gyrolab: Biotin & Alexafluor

SPARCL: HRP & Acridan



- Bridging ELISA
- Varying labels for capture and detection reagent
- Tocilizumab was US drug with varying lot numbers

MSD

GYROLAB

SPARCL

Prepare and Load Samples

Rexxip H

MRD = 1 in 10

Prepare Capture Solution

Capture: BioRad Clone, AbD21362 (Biotinylated)

350 nM, Bioaffy Bump Liquid

Prepare Detection Incubation


Detection: BioRad Clone, AbD21338_hlgG1 (AlexaFluor)

1.67 nM, Rexxip F

Analysis

Bioaffy 200 CD



Capture-Analyte-Detect (1% PMT)



MSD

GYROLAB

SPARCL





03

Data, Data, Data!



RANGE AND SENSITIVITY



- Gyrolab
- MSD
- ▲ SPARCL



Data normalised to signal to noise due to differing units.

- MSD: RLU
- Gyrolab: RFU
- SPARCL: RLU(AUC)




Range (LLOQ-ULOQ):

- MSD: 1.50 – 3,500 ng/mL
- Gyrolab: 15.0 – 20,000 ng/mL
- SPARCL: 400 – 7,500 ng/mL


INTER-ASSAY A&P

INTER-ASSAY A&P


Inter-Assay A&P - MSD



Inter-Assay A&P - SPARCL





Inter-Assay A&P - Gyrolab



- All assay met A&P criteria across their range
- Gyrolab demonstrated poor recovery at the ULOQ.
- MSD & SPARCL demonstrated increased TE in the low range.
- SPARCL had excellent recovery QCL-QCH

PROZONE/DILUTIONAL LINEARITY ASSESSMENT



- MSD
- Gyrolab
- ▲ SPARCL



Data normalised to signal to noise due to differing units.



MSD: RLU
Gyrolab: RFU
SPARCL: RLU(AUC)



Prozone Observed:

MSD: No
Gyrolab: No
SPARCL: Yes, >7500 ng/mL

PROZONE/DILUTIONAL LINEARITY ASSESSMENT



- MSD
- Gyrolab (Red = Failure)
- ▲ SPARCL



Data normalised to signal to noise due to differing units.

MSD: RLU

Gyrolab: RFU

SPARCL: RLU(AUC)




Dilutional Linearity was observed in all assays.

MSD: 1 in 400

Gyrolab: 1 in 1000

SPARCL: 1 in 1250

SELECTIVITY



Individual serum samples tested blank and at LLOQ.

No Individuals tested:

- MSD: 16
- Gyrolab: 10
- SPARCL: 14

Level	MSD	Gyrolab®	SPARCL™
Blank	100%	100%	100%
LLOQ	94%	100%	93%

04

Lessons learned and the future perspective.



Parameter	MSD	Gyrolab	SPARCL
Range (ng/mL)	1.50 – 3500	15.0 - 20,000	400 – 7500
MRD	1 in 10	1 in 10	1 in 10
A&P	Pass	Pass	Pass
Prozone	No	No	Extreme
Dilutional Linearity	1 in 400	1 in 1000	1 in 1250
Selectivity	Pass	Pass	Pass

TECHNOLOGICAL

OPERATIONAL


FUTURE

PERSPECTIVE

Parameter	MSD	Gyrolab	SPARCL
Run Time	5 Hours	2 Hours	1 Hour
No. Samples / Plate	32	40	32
No. Plates / Day	3	5*	12
Total Samples / Day**	96	200	384
*Utilising Gyrolab Xplore Model			
**Single Analyst			

OPERATIONAL

**FUTURE
PERSPECTIVE**



Future SPARCL PK work:

- Re-optimisation of master mix reagents and concentrations
- Re-optimisation of background reducing agent.

Future considerations:

- SPARCL technology and qualitative assays (ADA)
- Impact of automation on throughput.

**FUTURE
PERSPECTIVE**

Technology choice matters...

- MSD gives best sensitivity
- Gyrolab gives best dynamic range and least impact from matrix effects
- SPARCL gives the highest throughput but was restricted by prozone effect.

... but it must be matched with the critical reagents

- Critical reagent key quality attributes must be matched with the technology of interest to ensure the most appropriately designed assay.

MORGAN.EVANS@AGILEXBIOLABS.COM

Any questions?

