





TaMADOR
*Targeted Mass Spectrometry Assays
for Diabetes and Obesity Research*

An Antibody Cocktail–Based Immunoaffinity–LC–MS Method for Ultra-Sensitive and Robust Quantification of Circulating Proinsulin Forms and C-peptide

Wei-jun Qian and Jun Qu


Impaired proinsulin processing is a hallmark for β -cell dysfunction in diabetes



- ~70-90% loss of β -cell function at T1D diagnosis.
- ~40-50% loss at T2D diagnosis.

Diabetes 2021;70:1038–1050

Proinsulin forms and altered proinsulin processing



❑ Altered processing of proinsulin could be an important indicator of β -cell functions

- β -Cell Stress and Secretory Overload;
- Impaired Convertase Activity;
- ER Stress and Misfolding;
- Selective β -Cell Loss

➔ Selectively quantifying specific proinsulin forms in circulation is critical

Challenges in measuring proinsulin forms in circulation

Low levels in circulation

| Forms | Concentration in circulation |
|------------|---|
| Proinsulin | 10-200 pg/mL (healthy subjects) |
| Des-31,32 | Not previously quantified, but speculated to have comparable levels as proinsulin |
| Des-64,65 | Never detected in circulation before, might has much lower levels than Des-31,32 |
| C-peptide | 500-3000 pg/mL (healthy) |

Challenges in developing proinsulin form quantification assay

Severe cross-reactivity in most proinsulin ELISA kits used in clinical investigations


| Interference | Mercodia (Cat. 10-1118-01) | ALPCO (Cat. 80-PINHUT-CH01) |
|--------------|-------------------------------|--------------------------------|
| Des-31,32 | 95% | 100% |
| Des-64,65 | 84% | 100% |
| Insulin | <0.03% | <0.6% |
| C-peptide | <0.006% | <0.1% |

→ A method with high sensitivity and high selectivity is needed


Insufficient sensitivity in current LC-MS assays

| Reference | Sample preparation method | Instrumentation | Capability |
|-------------------------------------|---|---|---|
| J Chromatogr B, 2022, 1211, e123482 | <ul style="list-style-type: none"> • 100 μL plasma • PPT + SPE + online trap | <ul style="list-style-type: none"> • Micro-flow LC (3 μL/min) • Waters Xevo TQ-S | <ul style="list-style-type: none"> • Des-31,32: not quantitative; • C-peptide and insulin: quantitative; • Proinsulin: Undetectable in healthy subjects and some patients |
| Clin Chem 2021, 67, 854-862 | <ul style="list-style-type: none"> • 250 μL plasma • PPT + SPE | <ul style="list-style-type: none"> • High-flow LC • Thermo QE-Plus (high-resolution) | <ul style="list-style-type: none"> • Insulin and C-peptide: semi-quant • Proinsulin and des-31,32 were identified in patients • Undetectable proinsulin in healthy subjects • Des-64,65 was not detected |
| Forensic Toxicol, 2017, 35: 106-113 | <ul style="list-style-type: none"> • 100 μL plasma • PPT + immuno-enrichment + online trap | <ul style="list-style-type: none"> • Nano-flow LC (350 nL/min) • Thermo QE (high-resolution) | <ul style="list-style-type: none"> • LOD (LOQ) for proinsulin 50 pg/mL; • Proinsulin only detectable in 30% samples |

Efficient and uniform capture of proinsulin forms with an optimized antibody cocktail



Highly sensitive quantification with trapping-nano-LC-FAIMS/dCV-MS system



To produce selective surrogate peptides for specific forms

☐ Intact amino acid sequence

↓ Lys-C cleavage site



C-peptide:

EAEDLQVGQVELGGGPGAGSLQPLALEGSLQ



Proinsulin:

FVN...TPKTRREAEDLQVGQVELGGGPGAGSLQPLALEGSLQKRRGIV...YCN



Des-31,32:

EAEDLQVGQVELGGGPGAGSLQPLALEGSLQKRRGIV...YCN



Des-64,65:

FVN...TPKTRREAEDLQVGQVELGGGPGAGSLQPLALEGSLQ

To produce selective surrogate peptides for specific forms

- ❑ Lys-C generated surrogate peptides: idea by Wei-jun Qian



C-peptide: EAEDLQVGQVELGGPGAGSLQPLALEGSLQ

Proinsulin: TRREAEDLQVGQVELGGPGAGSLQPLALEGSLQK

Des-31,32: EAEDLQVGQVELGGPGAGSLQPLALEGSLQK

Des-64,65: TRREAEDLQVGQVELGGPGAGSLQPLALEGSLQ

Quantitative bias by pan-antibody capture across different forms

- Intact amino acid sequence



C-peptide:

EAEDLQVGQVELGGGPGAGSLQPLALEGSLQ

Proinsulin:

FVN...TPKTRREAEDLQVGQVELGGGPGAGSLQPLALEGSLQKRRGIV...YCN

Des-31,32:

EAEDLQVGQVELGGGPGAGSLQPLALEGSLQKRRGIV...YCN



Des-64,65:

FVN...TPKTRREAEDLQVGQVELGGGPGAGSLQPLALEGSLQ



Anti-C-peptide antibody (pan mAb)

Quantitative bias by pan-antibody capture across different forms



- PI:C ratio: molar ratio of proinsulin and C-peptide, multiplied by 100;
- Healthy PI:C=1~2; Diabetes PI:C= 2~10 in different publications.

Anal. Chem. 2025, in press;
doi.org/10.1021/acs.analchem.5c03256

Antibody-cocktail: improvement of capturing coverage

Patent Pending

B Chain

C-peptide

A Chain

C-peptide:

EAEDLQVGQVELGGGPGAGSLQPLALEGSLQ

Proinsulin:

FVN...TPKTRREAEDLQVGQVELGGGPGAGSLQPLALEGSLQKRRGIV...YCN

Des-31,32:

EAEDLQVGQVELGGGPGAGSLQPLALEGSLQKRRGIV...YCN

Des-64,65:

FVN...TPKTRREAEDLQVGQVELGGGPGAGSLQPLALEGSLQ




Uniform immunocapture with an antibody cocktail


(Clone C-pep-01+5E4/3+7F8)

A  **Single antibody**  **Antibody cocktail**


Single antibody **Antibody cocktail**



Spiked PI:C=0.1




Spiked PI:C=30



Anal. Chem. 2025, in press;
doi.org/10.1021/acs.analchem.5c03256

Trapping nano-LC FAIMS-dCV-MS system




Large ID trap:

- ✓ **Exceptional robustness;**
- ✓ High sensitivity;
- ✓ **Selective trap and delivery:**
remove matrix components;
- ✓ Rapid loading;
- ✓ Reproducible gradient delivery

Dynamic trapping/nano flows

- ✓ Rapid mass transfer without band broadening;
- ✓ High throughput: **cycle time~15 min;**
- ✓ High ESI response
- ✓ **Exceptional robustness**

FAIMS: separation mechanism orthogonal to LC or MS




Field Asymmetric Ion Mobility Spectrometry (FAIMS) VS conventional IMS

FAIMS: differential ion mobility



Conventional ion mobility spectrometry




Two advantages of FAIMS compared with conventional IMS

- **Continuous filter**, allows uninterrupted sampling of the ion stream (100% ion time utilization)
- **More orthogonal to m/z** , which can significantly improve the separation of analytes

FAIMS-dCV

Take full advantage of the extra selectivity by FAIMS in targeted protein quantification

Anal. Chem. 2024, 96, 32, 13140–13149




Peak CV:
Lower S/N



Optimal CV:
Higher S/N


A novel different CV (dCV) method for ultra-sensitive and selective protein quantification in plasma/tissues




B

Proinsulin **Des-31,32** **Des-64,65**

Without FAIMS




FAIMS/dCV




C

Exceptional robustness




No need for instrument cleaning/maintenance for >4 months


FAIMS significantly improves the intensities of peptide ions under nano-spray



De-clustering of adducted ions




Droplet evaporation




Calibration curves in surrogate matrix (equine serum)


C-peptide




Proinsulin



Des-31,32



Des-64,65



The range of concentration:

C-peptide: 8.8 – 4400 pg/mL;

Proinsulin: 1.7 – 340 pg/mL;

Des-31,32: 2.3 – 460 pg/mL;

Des-64,65: 3.6 – 360 pg/mL.

Validation

Validation results for quantification of the four targets in serum

| Analyte | Linearity range (pg/mL) | QC accuracy (%)/CV% (n=3) | | | |
|------------|----------------------------|---------------------------|---------------------|---------------------|----------------------|
| | | ^a LOQ | ^b Low-QC | ^b Mid-QC | ^b High-QC |
| C-peptide | 8.8- 4400 | 97.5/16.6 | 96.0/4.3 | 94.1/4.9 | 107.1/6.7 |
| Proinsulin | 1.7 – 337 | 100.2/9.6 | 97.9/4.2 | 88.8/1.0 | 105.1/3.2 |
| Des-31,32 | 2.3 – 468 | 93.3/9.9 | 98.7/13.8 | 106.3/9.6 | 108.4/6.9 |
| Des-64,65 | 3.6 – 365 | 87.0/1.2 | 106.7/11.6 | 105.6/14.3 | 98.5/3.6 |


^aThe validated LOQs for C-peptide, Proinsulin, Des-31,32 and Des-64,65 are respectively 8.8, 1.7, 2.3 and 3.6 pg/mL in serum.

^bThe levels of High, Mid, and Low QC in serum: C-peptide 3520.0, 440.0, and 22.0pg/mL; Proinsulin at 272.0, 34.0 and 4.3pg/mL; Des-64,65 at 292.0,36.5 and 4.6 pg/mL; Des-31,32 at 372.0, 46.5 and 5.8pg/mL

***Anal. Chem. 2025, in press;
doi.org/10.1021/acs.analchem.5c03256***

Application in measuring a clinical cohort of 80 subjects

The demographic characteristics of the four groups of clinical samples (N=20/group)





| | At-risk | Healthy Control#1 | New onset T1D | Healthy Control#2 |
|--|-----------|-------------------|---------------|-------------------|
| Gender: | 6/13 | 6/13 | 8/12 | 8/12 |
| Female/Male, N | | | | |
| Age (yr): mean (SD) | 13 (10.0) | 13 (9.8) | 12.4 (5.9) | 12.5 (6.0) |
| BMI (kg/m²): mean (SD) | 20 (6.8) | 20 (5.8) | 19.7 (4.7) | 20.2 (4.7) |

Application in measuring a clinical cohort of 80 subjects

Proinsulin and Des-31,32 were quantified in **all samples**;



Des-64,65 was quantified in **~70% samples**






Anal. Chem. 2025, in press;
doi.org/10.1021/acs.analchem.5c03256


LC-MS vs. ELISA (Merckodia 10-1118-01)




LC-MS vs. ELISA (TECO intact proinsulin assay)



Specific Proinsulin Form to C-peptide Ratio




Specific Proinsulin Form to C-peptide Ratio




2nd Application: Longitudinal Samples after Clinical T1D Diagnosis (Hybrid closed loop arm from the CLVer Trial)

Forlenza G.P., et al, JAMA, 2023.

Concentrations (pM)




- 25 patients followed at 13W, 26W, 39W, and 52w after T1D diagnosis.



Summary

- **The first method for accurate measurement of intact proinsulin, des-31,32 proinsulin, des-64,65 proinsulin, and C-peptide in circulation**
 - ✓ Achieved ultra-high sensitivity (2-5 pg/mL) with excellent robustness and reasonable throughput (15 min/cycle).
 - ✓ Enabled profiling of proinsulin processing *in vivo* with unprecedented detail
 - ✓ Showed superior sensitivity and selectivity to a total proinsulin ELISA method.
- **Novel observations:**
 - ✓ All forms of proinsulin and their ratios to C-peptide are abnormally elevated in the T1D setting, consistent with β -cell processing dysfunction.
 - ✓ The three forms of proinsulin showed unique patterns, which may reflect different types of β -cell stress and disease subtypes/outcomes.
 - ✓ Among the three proinsulin forms, the **des-31,32 proinsulin-to-C-peptide ratio** appeared to be the most robust discriminator of T1D.



Summary

- **The first method for accurate measurement of intact proinsulin, des-31,32 proinsulin, des-64,65 proinsulin, and C-peptide in circulation**
 - ✓ Achieved ultra-high sensitivity (2-5 pg/mL) with excellent robustness and reasonable throughput (15 min/cycle).
 - ✓ Enabled profiling of proinsulin processing *in vivo* with unprecedented detail
 - ✓ Showed superior sensitivity and selectivity to a total proinsulin ELISA method.
- **Novel observations:**
 - ✓ All forms of proinsulin and their ratios to C-peptide are abnormally elevated in the T1D setting, consistent with β -cell processing dysfunction.
 - ✓ The three forms of proinsulin showed unique patterns, which may reflect different types of β -cell stress and disease subtypes/outcomes.
 - ✓ Among the three proinsulin forms, the **des-31,32 proinsulin-to-C-peptide ratio** appeared to be the most robust discriminator of T1D.

Summary

- **Potential clinical application of the assay:**

- ✓ Staging, early detection and risk stratification;
- ✓ Refining the classification and phenotypes of the T1D;
- ✓ Monitor disease progress and therapeutic response;
- ✓ Characterization of β -cell dysfunction (e.g. detect subclinical β -cell dysfunction) and monitor β -cell-centric treatment effects.

- **Limitation of the assay**

- ✓ Signature peptide-based assay may not represent the exact molecular form of the targets;
- ✓ Unknown format or variants (e.g., truncation or modifications) containing the SP could interfere with the assay.

Acknowledgement

PNNL team:

Wei-jun Qian

TaiDu Lin

Lorenz Nierves

Shane Kelly

Samantha Powers

Tujin Shi

Thomas Fillmore

Panshak Dakup

& Others!



Qingqing Shen

Wang Cao

Huiting Tang

Fen Tu

Ming Zhang

Xiaoyu Zhu

Jie Pu

Chao Xue



TaMADOR

*Targeted Mass Spectrometry Assays
for Diabetes and Obesity Research*

UW school of Medicine

Cedars-Sinai Medical Center

Clinical Collaborators:

Indiana University:

Emily Sims

Carmella Evans-Molina

Gabriela S. F. Monaco

University of Minnesota

Antoinette Moran

University of British Columbia:

Bruce Verchere



National Institute of
Diabetes and Digestive
and Kidney Diseases